

Параметры для ввода в программу анализатора Beckman-Coulter AU 480

Parameters		Specific Test Parameters					
General	LIH	ISE	HbA1c		Calculated Test	Range	
Test Name: АпоВ ▾				Type: Serum ▾	Operation Yes ▾		
Sample Volume	<input type="text" value="1"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit			
Pre-Dilution Rate	<input type="text" value="1"/> ▾			Min.OD	<input type="text" value="-2.0"/>	Max.OD	<input type="text" value="2.5"/>
Rgt. Volume	R1(R1-1) <input type="text" value="100"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit			
				First Low	<input type="text" value="-2.0"/>	High	<input type="text" value="2.5"/>
				Last Low	<input type="text" value="-2.0"/>	High	<input type="text" value="2.5"/>
	R2(R2-1) <input type="text" value="20"/> μL	Dilution	<input type="text" value="0"/> μL	Dynamic Range Low	<input type="text" value="0"/>	High	<input type="text" value="2.5"/>
Wavelength	Pri <input type="text" value="340"/> nm ▾	Sec.	None ▾nm	Correlation Factor A	<input type="text" value="1"/>	B	<input type="text" value="0"/>
Method	END ▾			Factor for Maker A	<input type="text" value="1"/>	B	<input type="text" value="0"/>
Reaction Slope	<input type="text" value="+"/> ▾			Onboard Stability Period	<input type="text" value="999"/> Day	<input type="text" value=""/> Hour	
Measuring Point1 First	<input type="text" value="0"/>	Last	<input type="text" value="27"/>				
Measuring Point2 First	<input type="text" value="0"/>	Last	<input type="text" value="10"/>				
Linearity Limit	<input type="text" value=""/> %						
Lag Time Check	<input type="text" value=""/> ▾						

Parameters		Specific Test Parameters					
General	LIH	ISE	HbA1c		Calculated Test	Range	
Test Name: АпоВ ▾				Type: Serum ▾			
Value/Flag:	<input type="text" value=""/> ▾	Level L:	<input type="text" value=""/>	Level H:	<input type="text" value=""/>		
Specific Ranges:				Panic Value			
	Sex	Year	Month	Year	Month	Low	High
<input type="checkbox"/>	1. M ▾	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value="0.8"/>	<input type="text" value="1.55"/>
<input type="checkbox"/>	2. F ▾	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value="0.75"/>	<input type="text" value="1.5"/>
<input type="checkbox"/>	3. <input type="text" value=""/> ▾	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/>	4. <input type="text" value=""/> ▾	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/>	5. <input type="text" value=""/> ▾	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/>	6. <input type="text" value=""/> ▾	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
	7. No demographics					<input type="text" value="0.75"/>	<input type="text" value="1.55"/>
	8. Not within expected values					<input type="text" value="0.75"/>	<input type="text" value="1.55"/>
Unit	<input type="text" value="g/L"/>					Decimal Places	<input type="text" value="2"/>

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: АпоВ ▾				Type: Serum ▾	<input type="radio"/> Use Serum Cal.
Calibration Type:	6AB ▾	Formula:	Spline ▾	Counts:	<input type="text" value="2"/> ▾
<Calibrator Parameters>					
Calibrator	OD	Conc	Factor/OD range		Slope Check
			Low	High	None ▾
Point 1:	Saline	<input type="text" value="0"/>	<input type="text" value="-0.1"/>	<input type="text" value="2.5"/>	
Point 2:	АпоА1/В-1	*	<input type="text" value="-0.1"/>	<input type="text" value="2.5"/>	Allowance Range Check
Point 3:	АпоА1/В-2	*	<input type="text" value="-0.1"/>	<input type="text" value="2.5"/>	<input type="radio"/> Reagent Blank <input type="text" value=""/>
Point 4:	АпоА1/В-3	*	<input type="text" value="-0.1"/>	<input type="text" value="2.5"/>	<input type="radio"/> Calibration <input type="text" value=""/>
Point 5:	АпоА1/В-4	*	<input type="text" value="-0.1"/>	<input type="text" value="2.5"/>	Advanced Calibration
Point 6:	АпоА1/В-5	*	<input type="text" value="-0.1"/>	<input type="text" value="2.5"/>	Operation <input type="text" value=""/> ▾
Point 7:					Interval (RB/ACAL) <input type="text" value=""/> ▾
Point 8:					
Point 9:					
Point 10:					
<Point Cal. For		No. of Correction Points	<input type="text" value=""/> ▾	Use Master Curve	<input type="text" value=""/> ▾ <input type="radio"/> Lot Calibration
<Master Curve>					
Calibrator	OD	Conc	Low	High	Stability
Point 1:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	Reagent Blank <input type="text" value=""/> Day <input type="text" value=""/> Hour
Point 2:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	Calibration <input type="text" value=""/> Day <input type="text" value=""/> Hour
MB Type Factor:	<input type="text" value=""/>	1-Point Calibration Point	<input type="text" value=""/> ▾	<input type="radio"/> with Conc-0	

Диапазон нормальных значений указан в соответствии с рекомендациями производителя реагентов. При использовании единиц измерения, отличающихся от приведенных, убедитесь, что значения стандартов, контрольных материалов, диапазонов нормальных значений и линейности метода введены в этих же единицах.

*-вводится из паспорта к калибратору (TruCal Apo A1/B, 5 уровней).

Контроль по TruLab L (липидный) уровень 1 и уровень 2 или TruLab N и P.