

Альбумин в моче, (микроальбумин)

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Type: Photometric Version: 1
Number: Assay availability: Date: 18.02.2020
Run controls for onboard reagents by: Time: 11:01:08
Operator: ADMIN

Reaction definition Reagent / Sample Validity checks

Reaction mode:

Primary Secondary Read times
Wavelength: / Main: --
Last required read:
Absorbance range: -- Color correction: --
Sample blank type: Blank: --

Assay 13 of 54

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Type: Photometric Version: 1
Number: Assay availability: Date: 20.02.2020
Run controls for onboard reagents by: Time: 17:14:07
Operator: FSE

Reaction definition Reagent / Sample Validity checks

Reagent:
Diluent:
Diluent dispense mode:

R1 R2
Reagent volume:
Water volume:
Dispense mode:

Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution
<input type="text" value="Standart"/>	<input type="text" value="9.1"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	= 1:1.00	<input checked="" type="radio"/>
<input type="text" value="1:4"/>	<input type="text" value="20.0"/>	<input type="text" value="9.1"/>	<input type="text" value="80"/>	<input type="text"/>	= 1:4.00	<input type="radio"/>
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	= 2:1.00	<input type="radio"/>

Assay of

Done Cancel

?

Альбумин в моче, (микроальбумин)

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: MkAlb-u Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: TruCalMalb
Replicates: 2 [Range 1 - 3]

Calibrator level: Blank: Water Concentration: 0.0000
Cal 1: TruCalMalb1
Cal 2: TruCalMalb2
Cal 3: TruCalMalb3
Cal 4: TruCalMalb4
Cal 5: TruCalMalb5
Cal 6: None

Assay of

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: MkAlb-u Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator: TruCalMalb

Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	9.1			
Cal 1: TruCalMalb1	9.1			
Cal 2: TruCalMalb2	9.1			
Cal 3: TruCalMalb3	9.1			
Cal 4: TruCalMalb4	9.1			
Cal 5: TruCalMalb5	9.1			
Cal 6: None				

Assay of

Done Cancel

?

Альбумин в моче, (микроальбумин)

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: MкAlb-u Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibration intervals:

Full interval: 0 (hours)
Adjust interval: 0 (hours)

Calibration type:

Adjust type: 2 Point
Adjust level: Cal 1
Default ordering type: Full

Assay of

Done
Cancel

?

Единицы измерения - мг/л; знаков после запятой 2 (Configure result units).

Диапазон метода: 3 - 350 мг/л (Configure assay parameters\ Results).

Калибратор TruCal Albumine U/CSF, набор из 5 уровней, нулевой - вода (или физраствор).
Контроль по TruLab Urine уровень 1 и уровень 2.