

## Параметры для ввода в программу анализатора Beckman-Coulter AU 480

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <b>Cyst-C</b> ▾		Type: <b>Serum</b> ▾		Operation <b>Yes</b> ▾	
Sample Volume	<input type="text" value="1"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/> ▾			Min.OD	<input type="text" value="-2.0"/>
Rgt. Volume	R1(R1-1) <input type="text" value="90"/> μL	Dilution	<input type="text" value="0"/> μL	Max.OD	<input type="text" value="2.5"/>
				Reagent OD Limit	
				First Low	<input type="text" value="-2.0"/>
				High	<input type="text" value="2.5"/>
				Last Low	<input type="text" value="-2.0"/>
				High	<input type="text" value="2.5"/>
				Dynamic Range Low	<input type="text" value="0.1"/>
				High	<input type="text" value="8"/>
Wavelength	Pri <input type="text" value="520"/> nm ▾	Sec.	<input type="text" value="NONE"/> nm ▾	Correlation Factor A	<input type="text" value="1"/>
Method	<input type="text" value="FIXED"/> ▾			Factor for Maker A	<input type="text" value="1"/>
Reaction Slope	<input type="text" value="+"/> ▾			B	<input type="text" value="0"/>
Measuring Point1 First	<input type="text" value="12"/>	Last	<input type="text" value="24"/>	B	<input type="text" value="0"/>
Measuring Point2 First		Last			
Linearity Limit	<input type="text" value=""/> %			Onboard Stability Period	<input type="text" value="999"/> Day <input type="text" value=""/> Hour
Lag Time Check	<input type="text" value=""/> ▾				

  

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <b>Cyst-C</b> ▾		Type: <b>Serum</b> ▾			
Value/Flag:	<input type="text" value=""/> ▾	Level L:	<input type="text" value=""/>	Level H:	<input type="text" value=""/>
Specific Ranges:					
	From	To		Low	High
<input type="checkbox"/> 1.	Sex ▾	Year	Month	Year	Month
<input type="checkbox"/> 2.					
<input type="checkbox"/> 3.					
<input type="checkbox"/> 4.					
<input type="checkbox"/> 5.					
<input type="checkbox"/> 6.					
<input type="checkbox"/> 7.	No demographics			<input type="text" value="0.53"/>	<input type="text" value="1.02"/>
<input type="checkbox"/> 8.	Not within expected values			<input type="text" value="0.53"/>	<input type="text" value="1.02"/>
Unit	<input type="text" value="mg/l"/>		Decimal Places	<input type="text" value="2"/>	
Panic Value					
Low <input type="text" value=""/> High <input type="text" value=""/>					

  

Parameters		Calibration Parameters		
Calibrators	Calibration Specific	STAT Table Calibration		
General	ISE			
Test Name: <b>Cyst-C</b> ▾		Type: <b>Serum</b> ▾		<input type="radio"/> Use Serum Cal.
Calibration Type:	<input type="text" value="6AB"/> ▾	Formula:	<input type="text" value="Spline"/> ▾	Counts: <input type="text" value="2"/> ▾
<Calibrator Parameters>				
Calibrator	OD	Conc	Factor/OD range	Slope Check
			Low High	<input type="text" value="None"/> ▾
Point 1:	Saline	<input type="text" value="0"/>	<input type="text" value="-0.1"/> <input type="text" value="2.5"/>	
Point 2:	TruCal Cyst-C1	*	<input type="text" value="-0.1"/> <input type="text" value="2.5"/>	Allowance Range Check
Point 3:	TruCal Cyst-C2	*	<input type="text" value="-0.1"/> <input type="text" value="2.5"/>	<input type="radio"/> Reagent Blank <input type="text" value=""/>
Point 4:	TruCal Cyst-C3	*	<input type="text" value="-0.1"/> <input type="text" value="2.5"/>	<input type="radio"/> Calibration <input type="text" value=""/>
Point 5:	TruCal Cyst-C4	*	<input type="text" value="-0.1"/> <input type="text" value="2.5"/>	Advanced Calibration
Point 6:	TruCal Cyst-C5	*	<input type="text" value="-0.1"/> <input type="text" value="2.5"/>	Operation <input type="text" value=""/> ▾
Point 7:				Interval (RB/ACAL) <input type="text" value=""/> ▾
Point 8:				
Point 9:				
Point 10:				
<Point Cal. For Master Curve>				
	No. of Correction Points	<input type="text" value=""/> ▾	Use Master Curve	<input type="radio"/> Lot Calibration
Calibrator	OD	Conc	Low High	Stability
Point 1:				Reagent Blank <input type="text" value=""/> Day <input type="text" value=""/> Hour
Point 2:				Calibration <input type="text" value=""/> Day <input type="text" value=""/> Hour
MB Type Factor:	<input type="text" value=""/> 1-Point Calibration Point	<input type="text" value=""/> ▾	<input type="radio"/> with Conc-0	

Диапазон нормальных значений указан в соответствии с рекомендациями производителя.

При использовании единиц измерения, отличающихся от приведенных, убедитесь, что значения стандартов, контрольных материалов, диапазонов нормальных значений и линейности метода введены в этих же единицах.

\* -вводится из паспорта к калибратору (TruCal Cystatin C).

Контроль по TruLab Cystatin уровень 1 и уровень 2.