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Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <input type="text" value="C4"/>		Type: <input type="text" value="Serum"/>		Operation: <input type="text" value="Yes"/>	
Sample Volume	<input type="text" value="2.6"/>	μL	Dilution	<input type="text" value="0"/>	μL
Pre-Dilution Rate	<input type="text" value="1"/>		OD Limit	<input type="text" value="-2.0"/>	Max.OD <input type="text" value="2.5"/>
Rgt. Volume	R1(R1-1)	<input type="text" value="90"/>	μL	Dilution	<input type="text" value="0"/>
	R2(R2-1)	<input type="text" value="18"/>	μL	Dilution	<input type="text" value="0"/>
Wavelength	Pri	<input type="text" value="340"/>	nm	Sec.	<input type="text" value="800"/>
Method		<input type="text" value="END"/>			
Reaction Slope		<input type="text" value="+"/>			
Measuring Point1	First	<input type="text" value="0"/>	Last	<input type="text" value="27"/>	
Measuring Point2	First	<input type="text" value="0"/>	Last	<input type="text" value="10"/>	
Linearity Limit		<input type="text" value=""/>	%		
Lag Time Check		<input type="text" value=""/>			
Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <input type="text" value="C4"/>		Type: <input type="text" value="Serum"/>			
Value/Flag:	<input type="text" value=""/>	Level L:	<input type="text" value=""/>	Level H:	<input type="text" value=""/>
Specific Ranges:					
	From	To		Low	High
<input type="checkbox"/> 1.	Sex <input type="text" value=""/>	Year <input type="text" value=""/>	Month <input type="text" value=""/>	Year <input type="text" value=""/>	Month <input type="text" value=""/>
<input type="checkbox"/> 2.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 3.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 4.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 5.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 6.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 7.	No demographics			<input type="text" value="0.1"/>	<input type="text" value="0.4"/>
<input type="checkbox"/> 8.	Not within expected values			<input type="text" value="0.1"/>	<input type="text" value="0.4"/>
Unit	<input type="text" value="g/L"/>		Decimal Places	<input type="text" value="2"/>	
Panic Value					
		Low	<input type="text" value=""/>	High	<input type="text" value=""/>
Parameters		Calibration Parameters			
Calibrators	Calibration Specific		STAT Table Calibration		
General	ISE				
Test Name: <input type="text" value="C4"/>		Type: <input type="text" value="Serum"/>		Use Serum Cal. <input type="checkbox"/>	
Calibration Type: <input type="text" value="6AB"/>		Formula: <input type="text" value="Spline"/>		Counts: <input type="text" value="2"/>	
<Calibrator Parameters>					
Calibrator	OD	Conc	Factor/OD range	Low	High
Point 1:	TruCal Prot-0		*	-0.1	2.5
Point 2:	TruCal Prot-1		*	-0.1	2.5
Point 3:	TruCal Prot-2		*	-0.1	2.5
Point 4:	TruCal Prot-3		*	-0.1	2.5
Point 5:	TruCal Prot-4		*	-0.1	2.5
Point 6:	TruCal Prot-5		*	-0.1	2.5
Point 7:					
Point 8:					
Point 9:					
Point 10:					
Slope Check: <input type="text" value="None"/>					
Allowance Range Check					
		Reagent Blank		<input type="text" value=""/>	
		Calibration		<input type="text" value=""/>	
Advanced Calibration Operation: <input type="text" value=""/>					
Interval (RB/ACAL): <input type="text" value=""/>					
<Point Cal. For		No. of Correction Points	<input type="text" value=""/>	Use Master Curve	<input type="checkbox"/>
Master Curve>		OD Range		Lot Calibration <input type="checkbox"/>	
Calibrator	OD	Conc	Low	High	Stability
Point 1:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	Reagent Blank <input type="text" value=""/>
Point 2:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	Day <input type="text" value=""/>
MB Type Factor:		1-Point Calibration Point		<input type="checkbox"/>	
				with Conc-0 <input type="checkbox"/>	

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