

## Параметры для ввода в программу анализатора Olympus AU 400/640

Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:	Alb	<	>	Type:	Serum	Operation:	Yes			
Sample:	Volume	2.0	µL	Dilution	0	µL	Pre-Dilution Rate:	1		
Reagents:	R1 Volume	200	µL	Dilution	0	µL	Min OD	Max OD		
	R2 Volume	0	µL	Dilution	0	µL	Reagent OD limit:	L	H	
Wavelength:	Pri.	540	∇	Sec.	700	∇	First L	-0.1	First H	2.5
Method:	END		∇	Last L	-0.1	Last H	2.5			
Reaction slope:	+		∇	Dynamic Range:	L	2	H	60		
Measuring Point 1:	First	0	Last	27	Correlation Factor:	A	1	B	0	
Measuring Point 2:	First		Last		On-board stability period:	999				
Linearity:			%							
No Lag Time:			∇							

Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:	Alb	<	>	Type:	Serum					
Value/Flag:		∇	Level L:		Level H:					
Normal Ranges:	Age L	Age H		L	H					
	Sex	Year	Month	Year	Month					
o 1.										
o 2.										
o 3.										
o 4.										
o 5.										
o 6.										
7.	None Selected					35	52			
8.	Out of Range		L	H						
Panic Value:			L	H		Unit:	g/L			
						Decimal places:	1			

Calibration Specific										
General		ISE								
Test Name:	Alb	<	>	Type:	Serum					
Calibration Type:	AB	∇	Formula:	Y=AX+B	∇	Counts:	2	Process:	CONC	∇
Point 1:	Cal. No.	*	OD		CONC	**	Factor/OD-L	-99999	Factor/OD-H	99999
Point 2:										
Point 3:										
Point 4:										
Point 5:										
Point 6:										
Point 7:										
1-Point Cal. Point:		<input type="checkbox"/>	With CONC-0	Slope Check	None	∇	Advanced Calibration:		∇	
MB Type Factor:				Calibration Stability Period:	999					

Диапазон нормальных значений указан в соответствии с рекомендациями производителя.

\*\* -вводится из паспорта к калибратору

Калибровка линейная, АВ (калибратор TruCalU)

Контроль TruLabN и TruLabP.

\* -вводится в соответствии с настройками системы