

# . RF

ИммуноТD (Иммунотурбидиметрический)



## Параметры для ввода в программу анализатора Beckman-Coulter AU 480

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <input type="text" value="RF"/> < > Type: <input type="text" value="Serum"/> Operation <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="7"/> $\mu\text{L}$	Dilution	<input type="text" value="0"/> $\mu\text{L}$	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD	<input type="text" value="-2.0"/> Max.OD <input type="text" value="2.5"/>
Rgt. Volume	R1(R1-1) <input type="text" value="90"/> $\mu\text{L}$	Dilution	<input type="text" value="0"/> $\mu\text{L}$	Reagent OD Limit	
				First Low	<input type="text" value="-2.0"/> High <input type="text" value="2.5"/>
				Last Low	<input type="text" value="-2.0"/> High <input type="text" value="2.5"/>
R2(R2-1)	<input type="text" value="18"/> $\mu\text{L}$	Dilution	<input type="text" value="0"/> $\mu\text{L}$	Dynamic Range Low	<input type="text" value="1"/> High <input type="text" value="500"/>
Wavelength	Pri <input type="text" value="340"/> nm	Sec.	<input type="text" value="800"/> nm	Correlation Factor A	<input type="text" value="1"/> B <input type="text" value="0"/>
Method	<input type="text" value="END"/>			Factor for Maker A	<input type="text" value="1"/> B <input type="text" value="0"/>
Reaction Slope	<input type="text" value="+"/> %			Onboard Stability Period	<input type="text" value="999"/> Day <input type="text" value=""/> Hour
Measuring Point1	First <input type="text" value="0"/>	Last	<input type="text" value="27"/>		
Measuring Point2	First <input type="text" value="0"/>	Last	<input type="text" value="10"/>		
Linearity Limit	<input type="text" value=""/> %				
Lag Time Check	<input type="text" value=""/>				

  

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <input type="text" value="RF"/> < > Type: <input type="text" value="Serum"/>					
Value/Flag:	<input type="text" value=""/>	Level L:	<input type="text" value=""/>	Level H:	<input type="text" value=""/>
Specific Ranges:					
	From	To	Low	High	Panic Value
<input type="checkbox"/> 1.	Sex <input type="text" value=""/> Year <input type="text" value=""/> Month <input type="text" value=""/>	Year <input type="text" value=""/> Month <input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	Low <input type="text" value=""/> High <input type="text" value=""/>
<input type="checkbox"/> 2.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
<input type="checkbox"/> 3.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
<input type="checkbox"/> 4.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
<input type="checkbox"/> 5.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
<input type="checkbox"/> 6.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
<input type="checkbox"/> 7.	No demographics		<input type="text" value="0"/>	<input type="text" value="15"/>	
<input type="checkbox"/> 8.	Not within expected values		<input type="text" value="0"/>	<input type="text" value="15"/>	
Unit	<input type="text" value="ME/ml"/>		Decimal Places	<input type="text" value="1"/>	

  

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="RF"/> < > Type: <input type="text" value="Serum"/> Use Serum Cal. <input type="checkbox"/>					
Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	Counts:	<input type="text" value="2"/>
<Calibrator Parameters>					
Calibrator	OD	Conc	Factor/OD range	Slope Check	<input type="text" value="None"/>
			Low High		
Point 1:	Saline	0	-0.1 2.5		
Point 2:	TruCal RF -1	*	-0.1 2.5	Allowance Range Check	
Point 3:	TruCal RF -2	*	-0.1 2.5	<input type="checkbox"/> Reagent Blank <input type="text" value=""/>	
Point 4:	TruCal RF -3	*	-0.1 2.5	<input type="checkbox"/> Calibration <input type="text" value=""/>	
Point 5:	TruCal RF -4	*	-0.1 2.5	Advanced Calibration	
Point 6:	TruCal RF -5	*	-0.1 2.5	Operation	<input type="text" value=""/>
Point 7:				Interval (RB/ACAL)	<input type="text" value=""/>
Point 8:					
Point 9:					
Point 10:					
<Point Cal. For No. of Correction Points <input type="text" value=""/> Use Master Curve <input type="checkbox"/> Lot Calibration <input type="checkbox"/>					
Master Curve>					
Calibrator	OD	Conc	Low High	Stability	
Point 1:				Reagent Blank	<input type="text" value=""/> Day <input type="text" value=""/> Hour
Point 2:				Calibration	<input type="text" value=""/> Day <input type="text" value=""/> Hour
MB Type Factor:	<input type="text" value=""/>	1-Point Calibration Point <input type="checkbox"/>			<input type="checkbox"/> with Conc-0

Диапазон нормальных значений указан в соответствии с рекомендациями производителя.

При использовании единиц измерения, отличающихся от приведенных, убедитесь, что значения стандартов, контрольных материалов, диапазонов нормальных значений и линейности метода введены в этих же единицах.

\* -вводится из паспорта к калибратору (TruCal RF). Первая точка - физраствор.

Контроль по TruLab Protein уровень 1 и уровень 2.