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Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Type: Version:
Number: Assay availability: Date:
Run controls for onboard reagents by: Time:
Operator:

Reaction definition Reagent / Sample Validity checks

Reaction mode:

	Primary	Secondary	Read times	
Wavelength:	<input type="text" value="340"/>	<input type="text" value="548"/>	Main:	<input type="text" value="19"/> -- <input type="text" value="28"/>
Last required read:	<input type="text" value="30"/>		Flex:	<input type="text" value="18"/> -- <input type="text" value="23"/>
Absorbance range:	<input type="text" value="0.600"/> -- <input type="text" value="3.200"/>		Color correction:	<input type="text" value="13"/> -- <input type="text" value="16"/>
Sample blank type:	<input type="text" value="None"/>			

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Type: Version:
Number: Assay availability: Date:
Run controls for onboard reagents by: Time:
Operator:

Reaction definition Reagent / Sample Validity checks

Reagent:
Diluent:
Diluent dispense mode:

	R1	R2
Reagent volume:	<input type="text" value="160"/>	<input type="text" value="40"/>
Water volume:	<input type="text"/>	<input type="text"/>
Dispense mode:	<input type="text" value="Type 2"/>	<input type="text" value="Type 2"/>

Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution
<input type="text" value="Standart"/>	: <input type="text" value="16.0"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	= <input type="text" value="1:1.00"/>	<input checked="" type="radio"/>
<input type="text" value="1:10"/>	: <input type="text" value="1.6"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	= <input type="text" value="1:9.33"/>	<input type="radio"/>
<input type="text" value="1:20"/>	: <input type="text" value="10.0"/>	<input type="text" value="8.0"/>	<input type="text" value="90"/>	<input type="text"/>	= <input type="text" value="1:19.26"/>	<input type="radio"/>

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Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: AST Assay number: *** Date: 20.02.2020
Calibration method: Linear Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: TruCal_U Concentration: 0.0000

Blank: Water

Cal 1: TruCal_U1

Cal 2: None

Cal 3: None

Cal 4: None

Cal 5: None

Cal 6: None

Replicates: 2 [Range 1 - 3]

Assay of

Done

Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: AST Assay number: *** Date: 20.02.2020
Calibration method: Linear Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator:

Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	16.0			
Cal 1: TruCal_U1	16.0			
Cal 2: None				
Cal 3: None				
Cal 4: None				
Cal 5: None				
Cal 6: None				

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