

Параметры для ввода в программу анализатора Beckman-Coulter AU 680

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <input type="text" value="C4"/>		Type: <input type="text" value="Serum"/>		Operation <input type="text" value="Yes"/>	
Sample Volume	<input type="text" value="3.6"/>	μL	Dilution	<input type="text" value="0"/>	μL
Pre-Dilution Rate	<input type="text" value="1"/>		OD Limit	<input type="text" value="-2.0"/>	Min.OD
Rgt. Volume	<input type="text" value="125"/>	μL	Dilution	<input type="text" value="0"/>	μL
			Reagent OD Limit	<input type="text" value="2.5"/>	Max.OD
			First	<input type="text" value="-2.0"/>	Low
			Last	<input type="text" value="-2.0"/>	Low
			High	<input type="text" value="2.5"/>	High
			High	<input type="text" value="2.5"/>	High
			Dynamic Range Low	<input type="text" value="0"/>	Low
			High	<input type="text" value="0.9"/>	High
			Correlation Factor A	<input type="text" value="1"/>	A
			B	<input type="text" value="0"/>	B
			Factor for Maker	<input type="text" value="1"/>	A
			B	<input type="text" value="0"/>	B
Wavelength	<input type="text" value="340"/>	nm	Sec.	<input type="text" value="800"/>	nm
Method	<input type="text" value="END"/>				
Reaction Slope	<input type="text" value="+"/>		Onboard Stability Period	<input type="text" value="999"/>	Day
Measuring Point1 First	<input type="text" value="0"/>		Last	<input type="text" value="27"/>	Hour
Measuring Point2 First	<input type="text" value="0"/>		Last	<input type="text" value="10"/>	Hour
Linearity Limit	<input type="text" value=""/>	%			
Lag Time Check	<input type="text" value=""/>				

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <input type="text" value="C4"/>		Type: <input type="text" value="Serum"/>			
Value/Flag:	<input type="text" value=""/>	Level L:	<input type="text" value=""/>	Level H:	<input type="text" value=""/>
Specific Ranges:					
	From	To		Low	High
<input type="checkbox"/>	Sex	Year	Month	Year	Month
<input type="checkbox"/>	1.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/>	2.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/>	3.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/>	4.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/>	5.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/>	6.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/>	7. No demographics			<input type="text" value="0.1"/>	<input type="text" value="0.4"/>
<input type="checkbox"/>	8. Not within expected values			<input type="text" value="0.1"/>	<input type="text" value="0.4"/>
Unit	<input type="text" value="g/L"/>	Decimal Places	<input type="text" value="2"/>		
Panic Value					
Low <input type="text" value=""/> High <input type="text" value=""/>					

Parameters		Calibration Parameters		
Calibrators	Calibration Specific	STAT Table Calibration		
General	ISE			
Test Name: <input type="text" value="C4"/>		Type: <input type="text" value="Serum"/>		<input type="checkbox"/> Use Serum Cal.
Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	Counts: <input type="text" value="2"/>
<Calibrator Parameters>				
Calibrator	OD	Conc	Factor/OD range	Slope Check
			Low High	<input type="text" value="None"/>
Point 1:	Saline	0	-0.1 2.5	
Point 2:	TruCal Prot-1	*	-0.1 2.5	
Point 3:	TruCal Prot-2	*	-0.1 2.5	
Point 4:	TruCal Prot-3	*	-0.1 2.5	
Point 5:	TruCal Prot-4	*	-0.1 2.5	
Point 6:	TruCal Prot-5	*	-0.1 2.5	
Point 7:				
Point 8:				
Point 9:				
Point 10:				
Allowance Range Check				
<input type="checkbox"/> Reagent Blank				
<input type="checkbox"/> Calibration				
Advanced Calibration				
Operation <input type="text" value=""/>				
Interval (RB/ACAL) <input type="text" value=""/>				
<Point Cal. For No. of Correction Points <input type="text" value=""/> Use Master Curve <input type="checkbox"/> Lot Calibration <input type="checkbox"/>				
Master Curve				
Calibrator	OD	Conc	Low High	Stability
Point 1:				Reagent Blank <input type="text" value=""/> Day <input type="text" value=""/> Hour
Point 2:				Calibration <input type="text" value=""/> Day <input type="text" value=""/> Hour
MB Type Factor: <input type="text" value=""/> 1-Point Calibration Point <input type="checkbox"/> with Conc-0 <input type="checkbox"/>				

Диапазон нормальных значений указан в соответствии с рекомендациями производителя.

При использовании единиц измерения, отличающихся от приведенных, убедитесь, что значения стандартов, контрольных материалов, диапазонов нормальных значений и линейности метода введены в этих же единицах.

*-вводится в соответствии с настройками системы

** -вводится из паспорта к калибратору (TruCal Protein).

Контроль по TruLab Protein уровень 1 и уровень 2.