

Аполипопротеин А1

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect с8000 / с4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: ApoA1 Type: Photometric Version: 1
Number: *** Assay availability: Enabled Date: 18.02.2020
Run controls for onboard reagents by: Lot Time: 11:01:08
Operator: ADMIN

Reaction definition Reagent / Sample Validity checks

Reaction mode: End up

Primary Secondary Read times
Wavelength: 572 / None Main: 32 -- 33
Last required read: 33
Absorbance range: 0.000 -- 2.500 Color correction: --
Sample blank type: Self Blank: 14 -- 16

Assay 13 of 54

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: ApoA1 Type: Photometric Version: 1
Number: *** Assay availability: Enabled Date: 20.02.2020
Run controls for onboard reagents by: Lot Time: 17:14:07
Operator: FSE

Reaction definition Reagent / Sample Validity checks

Reagent: ApoA1
Diluent: Saline
Diluent dispense mode: Type 0

R1 R2
Reagent volume: 188 38
Water volume: --
Dispense mode: Type 2 Type 2

Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution
Standart	1.5				1:1.00	<input checked="" type="radio"/>
1:4	25.0	1.5	75		2:4.05	<input type="radio"/>
					2:1.00	<input type="radio"/>

Assay of

Done Cancel

?

Аполипопротеин А1

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: ApoA1 Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: TruCalApo Concentration: 0.0000

Blank: Water

Cal 1: TruCalApo1
Cal 2: TruCalApo2
Cal 3: TruCalApo3
Cal 4: TruCalApo4
Cal 5: TruCalApo5
Cal 6: None

Replicates: 2 [Range 1 - 3]

Assay of

Done

Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: ApoA1 Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator: TruCalApo

Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	1.5			
Cal 1: TruCalApo1	1.5			
Cal 2: TruCalApo2	1.5			
Cal 3: TruCalApo3	1.5			
Cal 4: TruCalApo4	1.5			
Cal 5: TruCalApo5	1.5			
Cal 6: None				

Assay of

Done

Cancel

?

Аполипопротеин А1

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: ApoA1 Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibration intervals:

Full interval: 0 (hours)
Adjust interval: 0 (hours)

Calibration type:

Adjust type: 2 Point
Adjust level: Cal 1
Default ordering type: Full

Assay of

Done
Cancel

?

Единицы измерения - г/л; знаков после запятой 3 (Configure result units)

Диапазон метода 0.002-2.5 г/мл (Configure assay parameters\Results).

Референсные интервалы: Муж.1.1-1.7 г/л, Жен.1.2-1.9 г/л (Configure assay parameters\Results).
При наличии ЛИС, не обязательно.

Диапазон нормальных значений указан в соответствии с рекомендациями производителя реагентов.

Калибратор TruCal ApoA/B, готовится 5 уровней, нулевой - вода (или физраствор).

Контроль по TruLab Lipid уровень 1 и уровень 2 или TruLab N и TruLab P.

Ниже приведен возможный вариант настройки авторазбавления калибратора на борту
Не проверено!

Аполипопротеин А1

ImmunoTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Architect с8000 / с4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: ApoA1 Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: TruCalApo Concentration: 0.0000

Blank: Water

Cal 1: TruCalApo1

Cal 2: None

Cal 3: None

Cal 4: None

Cal 5: None

Cal 6: None

Replicates: 2 [Range 1 - 3]

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: ApoA1 Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator: TruCalApo

Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	15	8.5	15	
Cal 1: TruCalApo1	3	8.5	102	
Cal 2: TruCalApo1	6	8.5	99	
Cal 3: TruCalApo1	9	8.5	96	
Cal 4: TruCalApo1	12	8.5	93	
Cal 5: TruCalApo1	15	8.5	90	
Cal 6: None				

Done Cancel

?