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Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <input type="text" value="Myogl"/>		Type: <input type="text" value="Serum"/>		Operation: <input type="text" value="Yes"/>	
Sample Volume	<input type="text" value="2.5"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD	<input type="text" value="-2.0"/>
Rgt. Volume	R1(R1-1) <input type="text" value="75"/> μL	Dilution	<input type="text" value="0"/> μL	Max.OD	<input type="text" value="2.5"/>
	R2(R2-1) <input type="text" value="25"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	
				First Low	<input type="text" value="-2.0"/>
				High	<input type="text" value="2.5"/>
				Last Low	<input type="text" value="-2.0"/>
				High	<input type="text" value="2.5"/>
Wavelength	Pri <input type="text" value="570"/> nm	Sec.	<input type="text" value="NONE"/>	Dynamic Range Low	<input type="text" value="5"/>
Method	<input type="text" value="FIXED"/>			High	<input type="text" value="600"/>
Reaction Slope	<input type="text" value="+"/> μm			Correlation Factor A	<input type="text" value="1"/>
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	Factor for Maker A	<input type="text" value="1"/>
Measuring Point2 First		Last		B	<input type="text" value="0"/>
Linearity Limit	<input type="text" value=""/>			B	<input type="text" value="0"/>
Lag Time Check	<input type="text" value=""/>			Onboard Stability Period	<input type="text" value="999"/> Day <input type="text" value=""/> Hour

Parameters		Specific Test Parameters				
General	LIH	ISE	HbA1c	Calculated Test	Range	
Test Name: <input type="text" value="Myogl"/>		Type: <input type="text" value="Serum"/>				
Value/Flag:	<input type="text" value=""/>	Level L:	<input type="text" value=""/>	Level H:	<input type="text" value=""/>	
Specific Ranges:						
	From	To		Low	High	
<input type="checkbox"/> 1.	Sex <input type="text" value=""/>	Year <input type="text" value=""/>	Month <input type="text" value=""/>	Year <input type="text" value=""/>	Month <input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 2.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 3.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 4.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 5.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 6.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 7.	No demographics			<input type="text" value="0"/>	<input type="text" value="70"/>	
<input type="checkbox"/> 8.	Not within expected values			<input type="text" value="0"/>	<input type="text" value="70"/>	
Unit	<input type="text" value="mkg/l"/>		Decimal Places	<input type="text" value="1"/>		

Parameters		Calibration Parameters		
Calibrators	Calibration Specific	STAT Table Calibration		
General	ISE			
Test Name: <input type="text" value="Myogl"/>		Type: <input type="text" value="Serum"/>		<input type="checkbox"/> Use Serum Cal.
Calibration Type:	<input type="text" value="5AB"/>	Formula:	<input type="text" value="Spline"/>	Counts: <input type="text" value="2"/>
<Calibrator Parameters>				
Calibrator	OD	Conc	Factor/OD range	Slope Check
			Low High	<input type="text" value="None"/>
Point 1:	Saline	0	-0.1 2.5	
Point 2:	TruCal Myogl 1	*	-0.1 2.5	Allowance Range Check
Point 3:	TruCal Myogl 2	*	-0.1 2.5	<input type="checkbox"/>
Point 4:	TruCal Myogl 3	*	-0.1 2.5	<input type="checkbox"/> Reagent Blank
Point 5:	TruCal Myogl 4	*	-0.1 2.5	<input type="checkbox"/> Calibration
Point 6:				
Point 7:				
Point 8:				Advanced Calibration
Point 9:				Operation
Point 10:				<input type="text" value=""/>
				Interval (RB/ACAL)
				<input type="text" value=""/>
<Point Cal. For	No. of Correction Points	<input type="text" value=""/>	Use Master Curve	<input type="checkbox"/>
Master Curve>			<input type="checkbox"/>	<input type="checkbox"/> Lot Calibration
Calibrator	OD	Conc	Low High	Stability
Point 1:				Reagent Blank
Point 2:				Day <input type="text" value=""/>
				Hour <input type="text" value=""/>
				Calibration
				Day <input type="text" value=""/>
				Hour <input type="text" value=""/>
MB Type Factor:	<input type="text" value=""/>	1-Point Calibration Point	<input type="checkbox"/>	<input type="checkbox"/> with Conc-0

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