

### Параметры для ввода в программу анализатора Beckman-Coulter AU 480

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <b>Bil-D</b> ▾		Type: <b>Serum</b> ▾		Operation <b>Yes</b> ▾	
Sample Volume	<input type="text" value="10"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/> ▾			Min.OD	<input type="text" value="-2.0"/>
Rgt. Volume	R1(R1-1) <input type="text" value="100"/> μL	Dilution	<input type="text" value="0"/> μL	Max.OD	<input type="text" value="2.5"/>
				Reagent OD Limit	
				First Low	<input type="text" value="-2.0"/>
				High	<input type="text" value="2.5"/>
				Last Low	<input type="text" value="-2.0"/>
				High	<input type="text" value="2.5"/>
R2(R2-1)	<input type="text" value="25"/> μL	Dilution	<input type="text" value="0"/> μL	Dynamic Range Low	<input type="text" value="0"/>
				High	<input type="text" value="170"/>
Wavelength	Pri <input type="text" value="540"/> nm ▾	Sec.	<input type="text" value="700"/> nm ▾	Correlation Factor A	<input type="text" value="1"/>
Method	<b>END</b> ▾			Factor for Maker A	<input type="text" value="1"/>
Reaction Slope	<input type="text" value="+"/> ▾			B	<input type="text" value="0"/>
Measuring Point1 First	<input type="text" value="0"/>	Last	<input type="text" value="27"/>	B	<input type="text" value="0"/>
Measuring Point2 First	<input type="text" value="0"/>	Last	<input type="text" value="10"/>	Onboard Stability Period	<input type="text" value="999"/> Day <input type="text" value=""/> Hour
Linearity Limit	<input type="text" value=""/> %				
Lag Time Check	<input type="text" value=""/> ▾				

  

Parameters		Specific Test Parameters				
General	LIH	ISE	HbA1c	Calculated Test	Range	
Test Name: <b>Bil-D</b> ▾		Type: <b>Serum</b> ▾				
Value/Flag:	<input type="text" value=""/> ▾	Level L:	<input type="text" value=""/>	Level H:	<input type="text" value=""/>	
Specific Ranges:						
	Sex	Year	Month	Year	Month	
<input type="checkbox"/>	1.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
<input type="checkbox"/>	2.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
<input type="checkbox"/>	3.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
<input type="checkbox"/>	4.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
<input type="checkbox"/>	5.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
<input type="checkbox"/>	6.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
	7. No demographics			<input type="text" value="0"/>	<input type="text" value="3.4"/>	
	8. Not within expected values			<input type="text" value="0"/>	<input type="text" value="3.4"/>	
Unit	<input type="text" value="mkmol/L"/>				Decimal Places	<input type="text" value="2"/>
Panic Value						
		Low	<input type="text" value=""/>	High	<input type="text" value=""/>	

  

Parameters		Calibration Parameters				
Calibrators	Calibration Specific	STAT Table Calibration				
General	ISE					
Test Name: <b>Bil-D</b> ▾		Type: <b>Serum</b> ▾		<input type="radio"/> Use Serum Cal.		
Calibration Type:	<b>AB</b> ▾	Formula:	<b>Y=AX+B</b> ▾	Counts:	<input type="text" value="2"/> ▾	
<Calibrator Parameters>						
Calibrator	OD	Conc	Factor/OD range	Slope Check	<input type="text" value="None"/> ▾	
Point 1:	<b>TruCal U</b>	<b>*</b>	<b>-99999 99999</b>	Allowance Range Check		
Point 2:				<input type="radio"/> Reagent Blank	<input type="text" value=""/>	
Point 3:				<input type="radio"/> Calibration	<input type="text" value=""/>	
Point 4:				Advanced Calibration		
Point 5:				Operation	<input type="text" value=""/> ▾	
Point 6:				Interval (RB/ACAL)	<input type="text" value=""/> ▾	
Point 7:				<Point Cal. For No. of Correction Points	<input type="text" value=""/> ▾	
Point 8:				Use Master Curve	<input type="text" value=""/> ▾	
Point 9:				<input type="radio"/> Lot Calibration		
Point 10:				Master Curve>		
Calibrator	OD	Conc	Low	High	Stability	
Point 1:					Reagent Blank	<input type="text" value=""/> Day <input type="text" value=""/> Hour
Point 2:					Calibration	<input type="text" value=""/> Day <input type="text" value=""/> Hour
MB Type Factor:	<input type="text" value=""/>	1-Point Calibration Point	<input type="text" value=""/> ▾	<input type="radio"/> with Conc-0		

Диапазон нормальных значений указан в соответствии с рекомендациями производителя. При использовании единиц измерения, отличающихся от приведенных, убедитесь, что значения стандартов, контрольных материалов, диапазонов нормальных значений и линейности метода введены в этих же единицах.

\*-вводится из паспорта к калибратору (калибратор TruCalU).

Контроль по TruLab N и TruLab P.