

. RF.

ИммуноТD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Beckman-Coulter AU 680

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <input type="text" value="RF"/>		Type: <input type="text" value="Serum"/>		Operation <input type="text" value="Yes"/>	
Sample Volume	<input type="text" value="9"/> $\mu\text{L}$	Dilution	<input type="text" value="0"/> $\mu\text{L}$	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD	<input type="text" value="-2.0"/> Max.OD <input type="text" value="2.5"/>
Rgt. Volume	R1(R1-1) <input type="text" value="125"/> $\mu\text{L}$	Dilution	<input type="text" value="0"/> $\mu\text{L}$	Reagent OD Limit	
				First Low	<input type="text" value="-2.0"/> High <input type="text" value="2.5"/>
				Last Low	<input type="text" value="-2.0"/> High <input type="text" value="2.5"/>
	R2(R2-1) <input type="text" value="25"/> $\mu\text{L}$	Dilution	<input type="text" value="0"/> $\mu\text{L}$	Dynamic Range Low <input type="text" value="1"/> High <input type="text" value="500"/>	
Wavelength	Pri <input type="text" value="340"/> $\text{nm}$	Sec.	<input type="text" value="800"/> $\text{nm}$	Correlation Factor A <input type="text" value="1"/> B <input type="text" value="0"/>	
Method	<input type="text" value="END"/>			Factor for Maker A <input type="text" value="1"/> B <input type="text" value="0"/>	
Reaction Slope	<input type="text" value="+"/>			Onboard Stability Period <input type="text" value="999"/> Day <input type="text" value=""/> Hour	
Measuring Point1	First <input type="text" value="0"/>	Last	<input type="text" value="27"/>		
Measuring Point2	First <input type="text" value="0"/>	Last	<input type="text" value="10"/>		
Linearity Limit	<input type="text" value=""/>	%			
Lag Time Check	<input type="text" value=""/>				

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <input type="text" value="RF"/>		Type: <input type="text" value="Serum"/>			
Value/Flag: <input type="text" value=""/>		Level L: <input type="text" value=""/>		Level H: <input type="text" value=""/>	
Specific Ranges:					
	From	To		Low	High
<input type="checkbox"/> 1.	Sex <input type="text" value=""/>	Year <input type="text" value=""/>	Month <input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 2.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 3.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 4.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 5.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 6.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
7. No demographics				<input type="text" value="0"/>	<input type="text" value="15"/>
8. Not within expected values				<input type="text" value="0"/>	<input type="text" value="15"/>
Unit <input type="text" value="ME/ml"/>		Decimal Places <input type="text" value="1"/>			

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="RF"/>		Type: <input type="text" value="Serum"/>		Use Serum Cal. <input type="checkbox"/>	
Calibration Type: <input type="text" value="6AB"/>		Formula: <input type="text" value="Spline"/>		Counts: <input type="text" value="2"/>	
<Calibrator Parameters>					
Calibrator	OD	Conc	Factor/OD range		Slope Check
			Low	High	<input type="text" value="None"/>
Point 1:	Saline	0	-0.1	2.5	
Point 2:	TruCal RF -1	*	-0.1	2.5	
Point 3:	TruCal RF -2	*	-0.1	2.5	
Point 4:	TruCal RF -3	*	-0.1	2.5	
Point 5:	TruCal RF -4	*	-0.1	2.5	
Point 6:	TruCal RF -5	*	-0.1	2.5	
Point 7:					
Point 8:					
Point 9:					
Point 10:					
<Point Cal. For		No. of Correction Points <input type="text" value=""/>	Use Master Curve <input type="checkbox"/>		Lot Calibration <input type="checkbox"/>
Master Curve>					
Calibrator	OD	Conc	Low	High	Stability
Point 1:					Reagent Blank <input type="text" value=""/> Day <input type="text" value=""/> Hour
Point 2:					Calibration <input type="text" value=""/> Day <input type="text" value=""/> Hour
MB Type Factor: <input type="text" value=""/>		1-Point Calibration Point <input type="checkbox"/>		with Conc-0 <input type="checkbox"/>	

Диапазон нормальных значений указан в соответствии с рекомендациями производителя.

При использовании единиц измерения, отличающихся от приведенных, убедитесь, что значения стандартов, контрольных материалов, диапазонов нормальных значений и линейности метода введены в этих же единицах.

\* -вводится из паспорта к калибратору (TruCal RF). Первая точка - физраствор.

Контроль по TruLab Protein уровень 1 и уровень 2.