

Параметры для ввода в программу анализатора Beckman-Coulter AU 680

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: Transf < >		Type: Serum < >		Operation: Yes < >	
Sample Volume	<input type="text" value="1.6"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/> < >			Min.OD	<input type="text" value="-2.0"/>
Rgt. Volume	R1(R1-1) <input type="text" value="200"/> μL	Dilution	<input type="text" value="0"/> μL	Max.OD	<input type="text" value="2.5"/>
				Reagent OD Limit	
				First Low	<input type="text" value="-2.0"/>
				High	<input type="text" value="2.5"/>
				Last Low	<input type="text" value="-2.0"/>
				High	<input type="text" value="2.5"/>
R2(R2-1)	<input type="text" value="40"/> μL	Dilution	<input type="text" value="0"/> μL	Dynamic Range Low	<input type="text" value="0"/>
				High	<input type="text" value="8"/>
Wavelength	Pri <input type="text" value="570"/> nm	Sec.	<input type="text" value="None"/> nm	Correlation Factor A	<input type="text" value="1"/>
Method	<input type="text" value="END"/> < >			Factor for Maker A	<input type="text" value="1"/>
Reaction Slope	<input type="text" value="+"/> < >			B	<input type="text" value="0"/>
Measuring Point1 First	<input type="text" value="0"/>	Last	<input type="text" value="27"/>	B	<input type="text" value="0"/>
Measuring Point2 First	<input type="text" value="0"/>	Last	<input type="text" value="10"/>		
Linearity Limit	<input type="text" value=""/> %			Onboard Stability Period	<input type="text" value="999"/> Day <input type="text" value=""/> Hour
Lag Time Check	<input type="text" value=""/> < >				

Parameters		Specific Test Parameters				
General	LIH	ISE	HbA1c	Calculated Test	Range	
Test Name: Transf < >		Type: Serum < >				
Value/Flag:	<input type="text" value=""/>	Level L:	<input type="text" value=""/>	Level H:	<input type="text" value=""/>	
Specific Ranges:						
	From	To		Low	High	
<input type="checkbox"/> 1.	Sex <input type="text" value=""/>	Year <input type="text" value=""/>	Month <input type="text" value=""/>	Year <input type="text" value=""/>	Month <input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 2.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 3.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 4.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 5.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 6.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 7.	No demographics			<input type="text" value="2"/>	<input type="text" value="3.6"/>	
<input type="checkbox"/> 8.	Not within expected values			<input type="text" value="2"/>	<input type="text" value="3.6"/>	
Unit	<input type="text" value="g/L"/>		Decimal Places	<input type="text" value="2"/>		

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: Transf < >		Type: Serum < >		<input type="radio"/> Use Serum Cal.	
Calibration Type:	<input type="text" value="6AB"/> < >	Formula:	<input type="text" value="Spline"/> < >	Counts:	<input type="text" value="2"/> < >
<Calibrator Parameters>					
Calibrator	OD	Conc	Factor/OD range	Slope Check	<input type="text" value="None"/> < >
Point 1:	Saline	<input type="text" value="0"/>	<input type="text" value="-0.1"/> <input type="text" value="2.5"/>		
Point 2:	TruCal Prot-1	<input type="text" value="*"/>	<input type="text" value="-0.1"/> <input type="text" value="2.5"/>	Allowance Range Check	
Point 3:	TruCal Prot-2	<input type="text" value="*"/>	<input type="text" value="-0.1"/> <input type="text" value="2.5"/>	<input type="radio"/> Reagent Blank	<input type="text" value=""/>
Point 4:	TruCal Prot-3	<input type="text" value="*"/>	<input type="text" value="-0.1"/> <input type="text" value="2.5"/>	<input type="radio"/> Calibration	<input type="text" value=""/>
Point 5:	TruCal Prot-4	<input type="text" value="*"/>	<input type="text" value="-0.1"/> <input type="text" value="2.5"/>	Advanced Calibration	
Point 6:	TruCal Prot-5	<input type="text" value="*"/>	<input type="text" value="-0.1"/> <input type="text" value="2.5"/>	Operation	<input type="text" value=""/> < >
Point 7:				Interval (RB/ACAL)	<input type="text" value=""/> < >
Point 8:					
Point 9:					
Point 10:					
<Point Cal. For		No. of Correction Points	<input type="text" value=""/> < >	Use Master Curve	<input type="text" value=""/> < >
Master Curve>		OD Range		<input type="radio"/> Lot Calibration	
Calibrator	OD	Conc	Low	High	Stability
Point 1:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	Reagent Blank
Point 2:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	Day
					Hour
					Day
					Hour
MB Type Factor:	<input type="text" value=""/>	1-Point Calibration Point	<input type="text" value=""/> < >	<input type="radio"/> with Conc-0	

Диапазон нормальных значений указан в соответствии с рекомендациями производителя.

При использовании единиц измерения, отличающихся от приведенных, убедитесь, что значения стандартов, контрольных материалов, диапазонов нормальных значений и линейности метода введены в этих же единицах.

* -вводится из паспорта к калибратору (TruCal Protein). Первая точка - физраствор.

Контроль по TruLab Protein уровень 1 и уровень 2 или по TruLab N и P.