

Параметры для ввода в программу анализатора iLab-Taurus

Test Name: NN Alb Test Code: IL_ Alb-d Measure Sample Reagent Ranges Limits Calibration

Sample Type: Serum

Reaction Cycle: Standard Extended

Reporting Unit: г/л * Decimal Points: 2 User Define Conv.Factor: 0.000

Methodology: Type: End Point Rate Measuring Point: 33 Photometric: 1 Wavelength 2 Wavelength Primary: 546 Secondary: 700

Correction Constant: Slope: 1.000 Intercept: 0.000

Test Name: NN Alb Test Code: IL_ Alb-d Measure Sample Reagent Ranges Limits Calibration

Sample Volume:

	Sample		Dilution	
	Volume	Sample Vol.	Diluent Vol.	
1	2.0	0.0	0.0	
2	2.0	20.0	80.0	
3	4.0	0.0	0.0	
4				

Diluent: Diluent Warning Limit: 0 Tests

Reagent Volume:

R1: Alb-d R2:

	Volume	Diluent Vol.	stiring	Warning Limit(tests)	Stability(days)
R1	200	0.0	<input checked="" type="checkbox"/>	20	No Control
R2			<input type="checkbox"/>		

Sampling Condition:

Condition No.	1	2	3	4
First Run	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Samp.Vol.Reduction	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Below N-Range	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Above N-Range	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Panic L	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Panic H	<input type="checkbox"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
User Range L	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/> User Range H	<input type="checkbox"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Noise	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Prozone	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/> HIGH!	<input type="checkbox"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/> ABS!	<input type="checkbox"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

Ranges			Limits																																																																																
<table border="1"> <thead> <tr> <th colspan="3">Normal Range</th> </tr> <tr> <th></th> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>Male</td> <td>35</td> <td>52</td> </tr> <tr> <td>Female</td> <td>35</td> <td>52</td> </tr> <tr> <td>Other</td> <td>35</td> <td>52</td> </tr> </tbody> </table>			Normal Range				Lower	Upper	Male	35	52	Female	35	52	Other	35	52	<table border="1"> <thead> <tr> <th>Reaction Slope</th> <th>Absorbance Limit</th> </tr> </thead> <tbody> <tr> <td> <input type="radio"/> Negative <input checked="" type="radio"/> Positive </td> <td> <input checked="" type="radio"/> Above <input type="radio"/> Below </td> </tr> <tr> <td> <table border="1"> <thead> <tr> <th>Non-Linear Limit</th> </tr> </thead> <tbody> <tr> <td>0 %</td> </tr> </tbody> </table> </td> <td> <table border="1"> <tbody> <tr> <td>2500.0 mAbs</td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th>Prozone Limit</th> </tr> </thead> <tbody> <tr> <td> <input type="radio"/> Above <input checked="" type="radio"/> Below </td> </tr> <tr> <td> <table border="1"> <thead> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th>Judge Point</th> </tr> </thead> <tbody> <tr> <td>0</td> </tr> </tbody> </table> </td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="3"> <table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>2</td> <td>60</td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> <tr> <td colspan="3"> <table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>2</td> <td>60</td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> <tr> <td colspan="3"> <p>Qualitative</p> <input type="radio"/> On <input checked="" type="radio"/> Off </td> <td colspan="2"></td> </tr> <tr> <td colspan="3"></td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table> </td> </tr> </tbody> </table>		Reaction Slope	Absorbance Limit	<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below	<table border="1"> <thead> <tr> <th>Non-Linear Limit</th> </tr> </thead> <tbody> <tr> <td>0 %</td> </tr> </tbody> </table>	Non-Linear Limit	0 %	<table border="1"> <tbody> <tr> <td>2500.0 mAbs</td> </tr> </tbody> </table>	2500.0 mAbs	<table border="1"> <thead> <tr> <th>Prozone Limit</th> </tr> </thead> <tbody> <tr> <td> <input type="radio"/> Above <input checked="" type="radio"/> Below </td> </tr> <tr> <td> <table border="1"> <thead> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th>Judge Point</th> </tr> </thead> <tbody> <tr> <td>0</td> </tr> </tbody> </table> </td> </tr> </tbody> </table>		Prozone Limit	<input type="radio"/> Above <input checked="" type="radio"/> Below	<table border="1"> <thead> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table>	Limit	Equation	0.0	none	<table border="1"> <thead> <tr> <th>Judge Point</th> </tr> </thead> <tbody> <tr> <td>0</td> </tr> </tbody> </table>		Judge Point	0	<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>2</td> <td>60</td> </tr> </tbody> </table>			User Range		Lower	Upper	2	60			<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>2</td> <td>60</td> </tr> </tbody> </table>			ValidRange		Lower	Upper	2	60			<p>Qualitative</p> <input type="radio"/> On <input checked="" type="radio"/> Off								<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>		Serum Index Limits			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0
Normal Range																																																																																			
	Lower	Upper																																																																																	
Male	35	52																																																																																	
Female	35	52																																																																																	
Other	35	52																																																																																	
Reaction Slope	Absorbance Limit																																																																																		
<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below																																																																																		
<table border="1"> <thead> <tr> <th>Non-Linear Limit</th> </tr> </thead> <tbody> <tr> <td>0 %</td> </tr> </tbody> </table>	Non-Linear Limit	0 %	<table border="1"> <tbody> <tr> <td>2500.0 mAbs</td> </tr> </tbody> </table>	2500.0 mAbs																																																																															
Non-Linear Limit																																																																																			
0 %																																																																																			
2500.0 mAbs																																																																																			
<table border="1"> <thead> <tr> <th>Prozone Limit</th> </tr> </thead> <tbody> <tr> <td> <input type="radio"/> Above <input checked="" type="radio"/> Below </td> </tr> <tr> <td> <table border="1"> <thead> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th>Judge Point</th> </tr> </thead> <tbody> <tr> <td>0</td> </tr> </tbody> </table> </td> </tr> </tbody> </table>		Prozone Limit	<input type="radio"/> Above <input checked="" type="radio"/> Below	<table border="1"> <thead> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table>	Limit	Equation	0.0	none	<table border="1"> <thead> <tr> <th>Judge Point</th> </tr> </thead> <tbody> <tr> <td>0</td> </tr> </tbody> </table>		Judge Point	0																																																																							
Prozone Limit																																																																																			
<input type="radio"/> Above <input checked="" type="radio"/> Below																																																																																			
<table border="1"> <thead> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table>	Limit	Equation	0.0	none																																																																															
Limit	Equation																																																																																		
0.0	none																																																																																		
<table border="1"> <thead> <tr> <th>Judge Point</th> </tr> </thead> <tbody> <tr> <td>0</td> </tr> </tbody> </table>		Judge Point	0																																																																																
Judge Point																																																																																			
0																																																																																			
<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>2</td> <td>60</td> </tr> </tbody> </table>			User Range		Lower	Upper	2	60																																																																											
User Range																																																																																			
Lower	Upper																																																																																		
2	60																																																																																		
<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>2</td> <td>60</td> </tr> </tbody> </table>			ValidRange		Lower	Upper	2	60																																																																											
ValidRange																																																																																			
Lower	Upper																																																																																		
2	60																																																																																		
<p>Qualitative</p> <input type="radio"/> On <input checked="" type="radio"/> Off																																																																																			
			<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>		Serum Index Limits			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0																																																																					
Serum Index Limits																																																																																			
	Limit																																																																																		
Hemolysis	0.0																																																																																		
Icterus	0.0																																																																																		
Lipemia	0.0																																																																																		

Test Name	Test Code	Measure	Sample Reagent	Ranges Limits	Calibration																																																
NN Alb	IL_ Alb-d																																																				
<table border="1"> <thead> <tr> <th colspan="2">Calibration</th> </tr> </thead> <tbody> <tr> <td>Method</td> <td>Curve Type</td> </tr> <tr> <td>1-point</td> <td>Linear</td> </tr> <tr> <td>Repeats</td> <td>Stability</td> </tr> <tr> <td>2</td> <td>days</td> </tr> <tr> <td>No</td> <td></td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Calibrator</th> </tr> <tr> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>TruCalU **</td> </tr> <tr> <td>2</td> <td></td> </tr> <tr> <td>3</td> <td></td> </tr> <tr> <td>4</td> <td></td> </tr> <tr> <td>5</td> <td></td> </tr> </tbody> </table> </td> </tr> </tbody> </table>		Calibration		Method	Curve Type	1-point	Linear	Repeats	Stability	2	days	No		<table border="1"> <thead> <tr> <th colspan="2">Calibrator</th> </tr> <tr> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>TruCalU **</td> </tr> <tr> <td>2</td> <td></td> </tr> <tr> <td>3</td> <td></td> </tr> <tr> <td>4</td> <td></td> </tr> <tr> <td>5</td> <td></td> </tr> </tbody> </table>		Calibrator		Calibrator	Conc	1	TruCalU **	2		3		4		5		<table border="1"> <thead> <tr> <th colspan="2">Condition</th> </tr> </thead> <tbody> <tr> <td><input checked="" type="checkbox"/></td> <td>Reagent Blank</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Auto Reagent Blank by Bottle</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Ask for calibration when reagent lot changes</td> </tr> </tbody> </table>		Condition		<input checked="" type="checkbox"/>	Reagent Blank	<input type="checkbox"/>	Auto Reagent Blank by Bottle	<input type="checkbox"/>	Ask for calibration when reagent lot changes	<table border="1"> <thead> <tr> <th colspan="2">Calibration/R-Blank Limit</th> </tr> </thead> <tbody> <tr> <td>R-Blank Limit</td> <td>2500.0 mAbs</td> </tr> <tr> <td>Cal Reps Range</td> <td>20.0 %</td> </tr> <tr> <td>Min Cal Reps</td> <td>200.0 mAbs</td> </tr> <tr> <td>Factor Change</td> <td>20 %</td> </tr> <tr> <td>M-Point Curve Fit</td> <td>0.0 %</td> </tr> </tbody> </table>		Calibration/R-Blank Limit		R-Blank Limit	2500.0 mAbs	Cal Reps Range	20.0 %	Min Cal Reps	200.0 mAbs	Factor Change	20 %	M-Point Curve Fit	0.0 %
Calibration																																																					
Method	Curve Type																																																				
1-point	Linear																																																				
Repeats	Stability																																																				
2	days																																																				
No																																																					
<table border="1"> <thead> <tr> <th colspan="2">Calibrator</th> </tr> <tr> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>TruCalU **</td> </tr> <tr> <td>2</td> <td></td> </tr> <tr> <td>3</td> <td></td> </tr> <tr> <td>4</td> <td></td> </tr> <tr> <td>5</td> <td></td> </tr> </tbody> </table>		Calibrator		Calibrator	Conc	1	TruCalU **	2		3		4		5																																							
Calibrator																																																					
Calibrator	Conc																																																				
1	TruCalU **																																																				
2																																																					
3																																																					
4																																																					
5																																																					
Condition																																																					
<input checked="" type="checkbox"/>	Reagent Blank																																																				
<input type="checkbox"/>	Auto Reagent Blank by Bottle																																																				
<input type="checkbox"/>	Ask for calibration when reagent lot changes																																																				
Calibration/R-Blank Limit																																																					
R-Blank Limit	2500.0 mAbs																																																				
Cal Reps Range	20.0 %																																																				
Min Cal Reps	200.0 mAbs																																																				
Factor Change	20 %																																																				
M-Point Curve Fit	0.0 %																																																				

Альбумин

Бромкрезоловый зеленый



* -Рекомендуемый параметр.

** - Вводится из паспорта калибратору.

Тип калибровки: Линейная (калибратор TruCalU) или по стандарту из набора.
Контроль по TruLab N и TruLab P.