

Параметры для ввода в программу анализатора Beckman-Coulter AU 480

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: Transf ▾		Type: Serum ▾		Operation Yes ▾	
Sample Volume	<input type="text" value="1"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/> ▾			Min.OD	<input type="text" value="-2.0"/>
Rgt. Volume	R1(R1-1) <input type="text" value="125"/> μL	Dilution	<input type="text" value="0"/> μL	Max.OD	<input type="text" value="2.5"/>
				Reagent OD Limit	
				First Low	<input type="text" value="-2.0"/>
				High	<input type="text" value="2.5"/>
				Last Low	<input type="text" value="-2.0"/>
				High	<input type="text" value="2.5"/>
				Dynamic Range Low	<input type="text" value="0"/>
				High	<input type="text" value="8"/>
Wavelength	Pri <input type="text" value="570"/> nm ▾	Sec.	<input type="text" value="None"/> nm ▾	Correlation Factor A	<input type="text" value="1"/>
Method	<input type="text" value="END"/> ▾			Factor for Maker A	<input type="text" value="1"/>
Reaction Slope	<input type="text" value="+"/> ▾			B	<input type="text" value="0"/>
Measuring Point1 First	<input type="text" value="0"/>	Last	<input type="text" value="27"/>	B	<input type="text" value="0"/>
Measuring Point2 First	<input type="text" value="0"/>	Last	<input type="text" value="10"/>		
Linearity Limit	<input type="text" value=""/> %			Onboard Stability Period	<input type="text" value="999"/> Day <input type="text" value=""/> Hour
Lag Time Check	<input type="text" value=""/> ▾				

Parameters		Specific Test Parameters					
General	LIH	ISE	HbA1c	Calculated Test	Range		
Test Name: Transf ▾		Type: Serum ▾					
Value/Flag:	<input type="text" value=""/> ▾	Level L:	<input type="text" value=""/>	Level H:	<input type="text" value=""/>		
Specific Ranges:							
	From	To		Low	High		
<input type="checkbox"/> 1.	Sex ▾	Year	Month	Year	Month	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 2.	<input type="text" value=""/> ▾	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 3.	<input type="text" value=""/> ▾	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 4.	<input type="text" value=""/> ▾	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 5.	<input type="text" value=""/> ▾	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 6.	<input type="text" value=""/> ▾	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
	7. No demographics				<input type="text" value="2"/>	<input type="text" value="3.6"/>	
	8. Not within expected values				<input type="text" value="2"/>	<input type="text" value="3.6"/>	
Unit	<input type="text" value="g/L"/>		Decimal Places	<input type="text" value="2"/>		Panic Value	
				Low	<input type="text" value=""/>	High	<input type="text" value=""/>

Parameters		Calibration Parameters				
Calibrators	Calibration Specific	STAT Table Calibration				
General	ISE					
Test Name: Transf ▾		Type: Serum ▾		Use Serum Cal. <input type="radio"/>		
Calibration Type:	<input type="text" value="6AB"/> ▾	Formula:	<input type="text" value="Spline"/> ▾	Counts:	<input type="text" value="2"/> ▾	
<Calibrator Parameters>						
Calibrator	OD	Conc	Low	High	Slope Check	
Point 1:	Saline	0	-0.1	2.5	<input type="text" value="None"/> ▾	
Point 2:	TruCal Prot-1	*	-0.1	2.5	Allowance Range Check	
Point 3:	TruCal Prot-2	*	-0.1	2.5		
Point 4:	TruCal Prot-3	*	-0.1	2.5		
Point 5:	TruCal Prot-4	*	-0.1	2.5		
Point 6:	TruCal Prot-5	*	-0.1	2.5		
Point 7:					<input type="radio"/> Reagent Blank <input type="text" value=""/>	
Point 8:					<input type="radio"/> Calibration <input type="text" value=""/>	
Point 9:					Advanced Calibration	
Point 10:					Operation <input type="text" value=""/> ▾	
					Interval (RB/ACAL)	<input type="text" value=""/> ▾
<Point Cal. For		No. of Correction Points	<input type="text" value=""/> ▾	Use Master Curve	<input type="text" value=""/> ▾	<input type="radio"/> Lot Calibration
<Master Curve>						
Calibrator	OD	Conc	Low	High	Stability	
Point 1:					Reagent Blank <input type="text" value=""/> Day <input type="text" value=""/> Hour	
Point 2:					Calibration <input type="text" value=""/> Day <input type="text" value=""/> Hour	
MB Type Factor:	<input type="text" value=""/>	1-Point Calibration Point	<input type="text" value=""/> ▾	<input type="radio"/> with Conc-0		

Диапазон нормальных значений указан в соответствии с рекомендациями производителя.

При использовании единиц измерения, отличающихся от приведенных, убедитесь, что значения стандартов, контрольных материалов, диапазонов нормальных значений и линейности метода введены в этих же единицах.

* -вводится из паспорта к калибратору (TruCal Protein). Первая точка - физраствор.

Контроль по TruLab Protein уровень 1 и уровень 2 или по TruLab N и P.