

Параметры для ввода в программу анализатора Beckman-Coulter AU 480

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <input type="text" value="Lp(a)"/> < > Type: <input type="text" value="Serum"/> Operation <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="1.5"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD	<input type="text" value="-2.0"/> Max.OD <input type="text" value="2.5"/>
Rgt. Volume	R1(R1-1) <input type="text" value="60"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	
				First Low	<input type="text" value="-2.0"/> High <input type="text" value="2.5"/>
				Last Low	<input type="text" value="-2.0"/> High <input type="text" value="2.5"/>
	R2(R2-1) <input type="text" value="30"/> μL	Dilution	<input type="text" value="0"/> μL	Dynamic Range Low	<input type="text" value="10"/> High <input type="text" value="1100"/>
Wavelength	Pri <input type="text" value="700"/> nm	Sec.	<input type="text" value="NONE"/> nm	Correlation Factor A	<input type="text" value="1"/> B <input type="text" value="0"/>
Method	<input type="text" value="FIXED"/>			Factor for Maker A	<input type="text" value="1"/> B <input type="text" value="0"/>
Reaction Slope	<input type="text" value="+"/>			Onboard Stability Period	<input type="text" value="999"/> Day <input type="text" value=""/> Hour
Measuring Point1	First <input type="text" value="13"/>	Last	<input type="text" value="27"/>		
Measuring Point2	First <input type="text" value=""/>	Last	<input type="text" value=""/>		
Linearity Limit	<input type="text" value=""/>				
Lag Time Check	<input type="text" value=""/>				

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <input type="text" value="Lp(a)"/> < > Type: <input type="text" value="Serum"/>					
Value/Flag:	<input type="text" value=""/>	Level L:	<input type="text" value=""/>	Level H:	<input type="text" value=""/>
Specific Ranges:					
	Sex	Year	Month	Year	Month
<input type="checkbox"/> 1.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 2.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 3.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 4.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 5.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 6.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 7.	No demographics			<input type="text" value="0"/>	<input type="text" value="300"/>
<input type="checkbox"/> 8.	Not within expected values			<input type="text" value=""/>	<input type="text" value=""/>
Unit	<input type="text" value="mg/L"/>		Decimal Places	<input type="text" value="1"/>	
Panic Value					
Low <input type="text" value=""/> High <input type="text" value=""/>					

Parameters		Calibration Parameters		
Calibrators	Calibration Specific	STAT Table Calibration		
General	ISE			
Test Name: <input type="text" value="Lp(a)"/> < > Type: <input type="text" value="Serum"/> Use Serum Cal. <input type="checkbox"/>				
Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	
<Calibrator Parameters>		Counts:	<input type="text" value="2"/>	
Calibrator	OD	Conc	Factor/OD range	Slope Check
			Low High	<input type="text" value=""/>
Point 1:	Saline	*	-0.1 2.5	Allowance Range Check
Point 2:	TruCal Lip(a)1	*	-0.1 2.5	
Point 3:	TruCal Lip(a)2	*	-0.1 2.5	
Point 4:	TruCal Lip(a)3	*	-0.1 2.5	
Point 5:	TruCal Lip(a)4	*	-0.1 2.5	
Point 6:	TruCal Lip(a)5	*	-0.1 2.5	
Point 7:				<input type="checkbox"/> Reagent Blank <input type="text" value=""/>
Point 8:				<input type="checkbox"/> Calibration <input type="text" value=""/>
Point 9:				Advanced Calibration
Point 10:				Operation <input type="text" value=""/>
				Interval (RB/ACAL) <input type="text" value=""/>
<Point Cal. For	No. of Correction Points	<input type="text" value=""/>	Use Master Curve	<input type="checkbox"/> Lot Calibration
Master Curve>				
Calibrator	OD	Conc	Low High	Stability
Point 1:				Reagent Blank <input type="text" value=""/> Day <input type="text" value=""/> Hour
Point 2:				Calibration <input type="text" value=""/> Day <input type="text" value=""/> Hour
MB Type Factor:	<input type="text" value=""/>	1-Point Calibration Point	<input type="checkbox"/>	<input type="checkbox"/> with Conc-0

Диапазон нормальных значений указан в соответствии с рекомендациями производителя.

При использовании единиц измерения, отличающихся от приведенных, убедитесь, что значения стандартов, контрольных материалов, диапазонов нормальных значений и линейности метода введены в этих же единицах.

*-вводится из паспорта к калибратору (TruCal Lp(a)).

Контроль по TruLab Lp(a) уровень 1 и уровень 2.