

## Параметры для ввода в программу анализатора Olympus AU 400/640

| Specific Test Parameters |           |     |     |          |     |       |                            |         |            |         |      |
|--------------------------|-----------|-----|-----|----------|-----|-------|----------------------------|---------|------------|---------|------|
| General                  |           | LIH | ISE | Range    |     |       |                            |         |            |         |      |
| Test Name:               | CK-MB ▾   |     |     | <        | >   | Type: | Serum ▾                    |         | Operation: | Yes ▾   |      |
| Sample:                  | Volume    | 6.0 | μL  | Dilution | 0   | μL    | Pre-Dilution Rate:         | 1       |            |         |      |
| Reagents:                | R1 Volume | 120 | μL  | Dilution | 0   | μL    | Min OD                     | Max OD  |            |         |      |
|                          | R2 Volume | 30  | μL  | Dilution | 0   | μL    | L                          | -0.1    | H          | 2.5     |      |
| Wavelength:              | Pri.      | 340 | ▾   | Sec.     | 540 | ▾     | Reagent OD limit:          | First L | -0.1       | First H | 2.5  |
| Method:                  | RATE ▾    |     |     |          |     |       | Last L                     | -0.1    | Last H     | 2.5     |      |
| Reaction slope:          | + ▾       |     |     |          |     |       | Dynamic Range:             | L       | 0          | H       | 2000 |
| Measuring Point 1:       | First     | 20  |     | Last     | 27  |       | Correlation Factor:        | A       | 1          | B       | 0    |
| Measuring Point 2:       | First     |     |     | Last     |     |       | On-board stability period: | 999     |            |         |      |
| Linearity:               | 15 %      |     |     |          |     |       |                            |         |            |         |      |
| No Lag Time:             | NO ▾      |     |     |          |     |       |                            |         |            |         |      |

| Specific Test Parameters |                  |       |       |          |       |       |          |   |                   |  |
|--------------------------|------------------|-------|-------|----------|-------|-------|----------|---|-------------------|--|
| General                  |                  | LIH   | ISE   | Range    |       |       |          |   |                   |  |
| Test Name:               | CK-MB ▾          |       |       | <        | >     | Type: | Serum ▾  |   |                   |  |
| Value/Flag:              |                  |       |       | Level L: |       |       | Level H: |   |                   |  |
| Normal Ranges:           |                  | Age L |       | Age H    |       | L     |          | H |                   |  |
|                          | Sex              | Year  | Month | Year     | Month |       |          |   |                   |  |
| <input type="radio"/>    | 1. Male ▾        |       |       |          |       |       |          |   |                   |  |
| <input type="radio"/>    | 2. Female ▾      |       |       |          |       |       |          |   |                   |  |
| <input type="radio"/>    | 3. ▾             |       |       |          |       |       |          |   |                   |  |
| <input type="radio"/>    | 4. ▾             |       |       |          |       |       |          |   |                   |  |
| <input type="radio"/>    | 5. ▾             |       |       |          |       |       |          |   |                   |  |
| <input type="radio"/>    | 6. ▾             |       |       |          |       |       |          |   |                   |  |
|                          | 7. None Selected |       |       |          |       |       | 0        |   | 24                |  |
|                          | 8. Out of Range  |       |       |          |       |       |          |   |                   |  |
| Panic Value:             | L                |       | H     |          | Unit: |       | U/L      |   | Decimal places: 1 |  |

| Calibration Specific |          |                          |             |             |          |       |                               |     |             |          |        |
|----------------------|----------|--------------------------|-------------|-------------|----------|-------|-------------------------------|-----|-------------|----------|--------|
| General              |          | ISE                      |             |             |          |       |                               |     |             |          |        |
| Test Name:           | CK-MB ▾  |                          |             | <           | >        | Type: | Serum ▾                       |     |             |          |        |
| Calibration Type:    | MB ▾     |                          |             | Formula:    | Y=AX+B ▾ |       | Counts:                       |     |             | Process: | CONC ▾ |
| Point 1:             | Cal. No. |                          | OD          |             | CONC     |       | Factor/OD-L                   |     | Factor/OD-H |          |        |
| Point 2:             |          |                          |             |             |          |       |                               |     |             |          |        |
| Point 3:             |          |                          |             |             |          |       |                               |     |             |          |        |
| Point 4:             |          |                          |             |             |          |       |                               |     |             |          |        |
| Point 5:             |          |                          |             |             |          |       |                               |     |             |          |        |
| Point 6:             |          |                          |             |             |          |       |                               |     |             |          |        |
| Point 7:             |          |                          |             |             |          |       |                               |     |             |          |        |
| 1-Point Cal. Point:  |          | <input type="checkbox"/> | With CONC-0 | Slope Check | None ▾   |       | Advanced Calibration:         |     |             | ▾        |        |
| MB Type Factor:      | 8300     |                          |             |             |          |       | Calibration Stability Period: | 999 |             |          |        |

Диапазон нормальных значений указан в соответствии с рекомендациями производителя реагентов к набору и может быть изменен пользователем в соответствии с местными требованиями. Калибровка линейная, АВ (калибратор TruCal CK-MB <на регистрации>) или по фактору -МВ (фактор см. инструкцию; вводится в поле MB Type Factor). Контроль TruLab N и TruLab P.

Калибровка по фактору.