

## Параметры для ввода в программу анализатора iLab-Taurus

Test Name: NN UA Test Code: IL\_ UA-d

Measure Sample Reagent Ranges Limits Calibration

Sample Type: Serum

Reaction Cycle: Standard Extended

Reporting Unit: МКОЛЬ/л \* Decimal Points: 1

User Define: Conv.Factor: 0.000

Methodology: Type: End Point Rate

Measuring Point: 16 - 33

Photometric: 1 Wavelength 2 Wavelength

Primary: 546 Secondary: 700

Correction Constant: Slope: 1.000 Intercept: 0.000

Test Name: NN UA Test Code: IL\_ UA-d

Measure Sample Reagent Ranges Limits Calibration

Sample Volume: Sampling Condition:

Sample	Dilution	Diluent	Condition No.	1	2	3	4
1	3.25	0.0	First Run				
2	2.0	0.0	Samp.Vol.Reduction				
3	6.5	0.0	Below N-Range				
4	0.0	0.0	Above N-Range				
			Panic L				
			Panic H				
			User Range L				
			User Range H				
			Noise				
			Prozone				
			HIGH!				
			ABS!				

Diluent Warning Limit: 0 Tests

Reagent Volume: R1: UA-d1 R2: UA-d2

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	160	0.0		20	No Control
R2	40	0.0		20	No Control

# Мочевая кислота

TOOS (метод с толуидином (TOOS))

Ranges		Limits																																		
<table border="1"><thead><tr><th colspan="3">Normal Range</th></tr><tr><th></th><th>Lower</th><th>Upper</th></tr></thead><tbody><tr><td>Male</td><td>208</td><td>428</td></tr><tr><td>Female</td><td>155</td><td>357</td></tr><tr><td>Other</td><td>155</td><td>428</td></tr></tbody></table>		Normal Range				Lower	Upper	Male	208	428	Female	155	357	Other	155	428	<table border="1"><thead><tr><th>Reaction Slope</th><th>Absorbance Limit</th></tr></thead><tbody><tr><td><input type="radio"/> Negative <input checked="" type="radio"/> Positive</td><td><input checked="" type="radio"/> Above <input type="radio"/> Below</td></tr><tr><td><table border="1"><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>0</td><td>%</td></tr></thead><tbody></tbody></table></td><td><table border="1"><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td>Limit: 0.0</td><td>Equation: none</td></tr><tr><td colspan="2">Judge Point: 0</td></tr></thead></table></td></tr></tbody></table>		Reaction Slope	Absorbance Limit	<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below	<table border="1"><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>0</td><td>%</td></tr></thead><tbody></tbody></table>	Non-Linear Limit		0	%	<table border="1"><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td>Limit: 0.0</td><td>Equation: none</td></tr><tr><td colspan="2">Judge Point: 0</td></tr></thead></table>	Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below		Limit: 0.0	Equation: none	Judge Point: 0	
Normal Range																																				
	Lower	Upper																																		
Male	208	428																																		
Female	155	357																																		
Other	155	428																																		
Reaction Slope	Absorbance Limit																																			
<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below																																			
<table border="1"><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>0</td><td>%</td></tr></thead><tbody></tbody></table>	Non-Linear Limit		0	%	<table border="1"><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td>Limit: 0.0</td><td>Equation: none</td></tr><tr><td colspan="2">Judge Point: 0</td></tr></thead></table>	Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below		Limit: 0.0	Equation: none	Judge Point: 0																								
Non-Linear Limit																																				
0	%																																			
Prozone Limit																																				
<input type="radio"/> Above <input checked="" type="radio"/> Below																																				
Limit: 0.0	Equation: none																																			
Judge Point: 0																																				

	User Range			------------	-------		Lower	Upper		18	2500				
	ValidRange			------------	-------		Lower	Upper		18	2500				
Qualitative															
☐ On ☒ Off															

Test Name		Test Code		Measure		Sample Reagent		Ranges Limits		Calibration																								
NN	UA	IL	UA-d																															
Calibration		Condition				Calibration/R-Blank Limit																												
<table border="1"><thead><tr><th colspan="2">Method</th></tr><tr><td>1-point</td><td></td></tr></thead><tbody></tbody></table>		Method		1-point		<table border="1"><thead><tr><th colspan="2">Curve Type</th></tr><tr><td>Linear</td><td></td></tr></thead><tbody></tbody></table>		Curve Type		Linear		<input checked="" type="checkbox"/> Reagent Blank				<table border="1"><thead><tr><th colspan="2">R-Blank Limit</th></tr><tr><td>2500.0</td><td>mAbs</td></tr></thead><tbody></tbody></table>				R-Blank Limit		2500.0	mAbs											
Method																																		
1-point																																		
Curve Type																																		
Linear																																		
R-Blank Limit																																		
2500.0	mAbs																																	
<table border="1"><thead><tr><th colspan="2">Repeats</th></tr><tr><td>2</td><td></td></tr></thead><tbody></tbody></table>		Repeats		2		<table border="1"><thead><tr><th colspan="2">Stability</th></tr><tr><td></td><td>days</td></tr></thead><tbody></tbody></table>		Stability			days	<input type="checkbox"/> Auto Reagent Blank by Bottle				<table border="1"><thead><tr><th colspan="2">Cal Reps Range</th></tr><tr><td>20.0</td><td>%</td></tr></thead><tbody></tbody></table>				Cal Reps Range		20.0	%											
Repeats																																		
2																																		
Stability																																		
	days																																	
Cal Reps Range																																		
20.0	%																																	
<table border="1"><thead><tr><th colspan="2">No</th></tr><tr><td></td><td></td></tr></thead><tbody></tbody></table>		No				<input type="checkbox"/> Ask for calibration when reagent lot changes				<table border="1"><thead><tr><th colspan="2">Min Cal Reps</th></tr><tr><td>50.0</td><td>mAbs</td></tr></thead><tbody></tbody></table>				Min Cal Reps		50.0	mAbs																	
No																																		
Min Cal Reps																																		
50.0	mAbs																																	
<table border="1"><thead><tr><th colspan="3">Calibrator</th></tr><tr><th></th><th>Calibrator</th><th>Conc</th></tr></thead><tbody><tr><td>1</td><td>TruCalU</td><td>**</td></tr><tr><td>2</td><td></td><td></td></tr><tr><td>3</td><td></td><td></td></tr><tr><td>4</td><td></td><td></td></tr><tr><td>5</td><td></td><td></td></tr></tbody></table>		Calibrator				Calibrator	Conc	1	TruCalU	**	2			3			4			5							<table border="1"><thead><tr><th colspan="2">Factor Change</th></tr><tr><td>20</td><td>%</td></tr></thead><tbody></tbody></table>				Factor Change		20	%
Calibrator																																		
	Calibrator	Conc																																
1	TruCalU	**																																
2																																		
3																																		
4																																		
5																																		
Factor Change																																		
20	%																																	
						<table border="1"><thead><tr><th colspan="2">M-Point Curve Fit</th></tr><tr><td>0.0</td><td>%</td></tr></thead><tbody></tbody></table>				M-Point Curve Fit		0.0	%																					
M-Point Curve Fit																																		
0.0	%																																	

# Мочевая кислота

TOOS (метод с толуидином (TOOS))



\* -Рекомендуемый параметр.

\*\* - Вводится из паспорта калибратору.

Тип калибровки: Линейная (калибратор TruCalU) или по стандарту из набора.  
Контроль по TruLab N и TruLab P.