

# С-Реактивный белок, ЦРБ (CRP-FS)

ИммуноТД (Иммунотурбидиметрический)

## Параметры для ввода в программу анализатора iLab-Taurus

Test Name: NN **CRP** Test Code: IL\_ **CRP-d** **Measure** **Sample Reagent** **Ranges Limits** **Calibration**

Sample Type: **Serum**

Reaction Cycle: ☒ Standard ☐ Extended

Reporting Unit: ☒ Мг/л\* ☐ User Define   
Decimal Points: 2   
Conv. Factor: 0.000

Methodology:   
Type: ☒ End Point ☐ Rate   
Measuring Point: 16 - ☒ 33   
Photometric: ☐ 1 Wavelength ☒ 2 Wavelength   
Primary: 340 Secondary: 700   
Correction Constant:   
Slope: 1.000 Intercept: 0.000

Test Name: NN **CRP** Test Code: IL\_ **CRP-d** **Measure** **Sample Reagent** **Ranges Limits** **Calibration**

Sample Volume: 

	Sample	Dilution	
	Volume	Sample Vol.	Diluent Vol.
1	9.0	0.0	0.0
2	2.0	0.0	0.0
3	12.0	0.0	0.0
4	0.0	0.0	0.0

  
Diluent:   
Diluent Warning Limit: 0 Tests

Reagent Volume:   
R1: **CRP-1** R2: **CRP-2**   

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	150	0.0	<input checked="" type="checkbox"/>	20	No Control
R2	30	0.0	<input checked="" type="checkbox"/>	20	No Control

Sampling Condition: 

	Condition No.	1	2	3	4
<input type="checkbox"/>	First Run	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Samp.Vol.Reduction	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Below N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Above N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Panic L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Panic H	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	User Range L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	User Range H	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Noise	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Prozone	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	HIGH!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	ABS!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

# С-Реактивный белок, ЦРБ (CRP-FS)

ИммуноТД (Иммунотурбидиметрический)

Ranges	Limits																																													
<table><thead><tr><th colspan="3">Normal Range</th></tr><tr><th></th><th>Lower</th><th>Upper</th></tr></thead><tbody><tr><td>Male</td><td>0.0</td><td>0.5</td></tr><tr><td>Female</td><td>0.0</td><td>0.5</td></tr><tr><td>Other</td><td>0.0</td><td>0.5</td></tr></tbody></table>	Normal Range				Lower	Upper	Male	0.0	0.5	Female	0.0	0.5	Other	0.0	0.5	<table><thead><tr><th>Reaction Slope</th><th>Absorbance Limit</th></tr></thead><tbody><tr><td><input type="radio"/> Negative <input checked="" type="radio"/> Positive</td><td><input checked="" type="radio"/> Above <input type="radio"/> Below</td></tr><tr><td><table><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>100</td><td>%</td></tr></thead></table></td><td><table><tr><td>2600.0</td><td>mAbs</td></tr></table></td></tr><tr><td></td><td><table><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td><table><tr><td>Limit</td><td>Equation</td></tr><tr><td>0.0</td><td>none</td></tr></table></td><td></td></tr><tr><td><table><thead><tr><th colspan="2">Judge Point</th></tr><tr><td>0</td><td></td></tr></thead></table></td><td></td></tr></thead></table></td></tr></tbody></table>	Reaction Slope	Absorbance Limit	<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below	<table><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>100</td><td>%</td></tr></thead></table>	Non-Linear Limit		100	%	<table><tr><td>2600.0</td><td>mAbs</td></tr></table>	2600.0	mAbs		<table><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td><table><tr><td>Limit</td><td>Equation</td></tr><tr><td>0.0</td><td>none</td></tr></table></td><td></td></tr><tr><td><table><thead><tr><th colspan="2">Judge Point</th></tr><tr><td>0</td><td></td></tr></thead></table></td><td></td></tr></thead></table>	Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below		<table><tr><td>Limit</td><td>Equation</td></tr><tr><td>0.0</td><td>none</td></tr></table>	Limit	Equation	0.0	none		<table><thead><tr><th colspan="2">Judge Point</th></tr><tr><td>0</td><td></td></tr></thead></table>	Judge Point		0		
Normal Range																																														
	Lower	Upper																																												
Male	0.0	0.5																																												
Female	0.0	0.5																																												
Other	0.0	0.5																																												
Reaction Slope	Absorbance Limit																																													
<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below																																													
<table><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>100</td><td>%</td></tr></thead></table>	Non-Linear Limit		100	%	<table><tr><td>2600.0</td><td>mAbs</td></tr></table>	2600.0	mAbs																																							
Non-Linear Limit																																														
100	%																																													
2600.0	mAbs																																													
	<table><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td><table><tr><td>Limit</td><td>Equation</td></tr><tr><td>0.0</td><td>none</td></tr></table></td><td></td></tr><tr><td><table><thead><tr><th colspan="2">Judge Point</th></tr><tr><td>0</td><td></td></tr></thead></table></td><td></td></tr></thead></table>	Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below		<table><tr><td>Limit</td><td>Equation</td></tr><tr><td>0.0</td><td>none</td></tr></table>	Limit	Equation	0.0	none		<table><thead><tr><th colspan="2">Judge Point</th></tr><tr><td>0</td><td></td></tr></thead></table>	Judge Point		0																															
Prozone Limit																																														
<input type="radio"/> Above <input checked="" type="radio"/> Below																																														
<table><tr><td>Limit</td><td>Equation</td></tr><tr><td>0.0</td><td>none</td></tr></table>	Limit	Equation	0.0	none																																										
Limit	Equation																																													
0.0	none																																													
<table><thead><tr><th colspan="2">Judge Point</th></tr><tr><td>0</td><td></td></tr></thead></table>	Judge Point		0																																											
Judge Point																																														
0																																														

☐ On ☒ Off

Test Name	Test Code	Measure	Sample Reagent	Ranges Limits	Calibration
NN	CRP	IL_	CRP-d		

Calibration	Condition	Calibration/R-Blank Limit																									
<table><tr><td>Method</td><td>Curve Type</td></tr><tr><td>M-point</td><td>5P Log</td></tr><tr><td>Repeats</td><td>Stability</td></tr><tr><td>2</td><td></td></tr><tr><td>No</td><td></td></tr><tr><td>6</td><td></td></tr></table>	Method	Curve Type	M-point	5P Log	Repeats	Stability	2		No		6		<table><tr><td><input type="checkbox"/> Reagent Blank</td></tr><tr><td><input type="checkbox"/> Auto Reagent Blank by Bottle</td></tr><tr><td><input type="checkbox"/> Ask for calibration when reagent lot changes</td></tr></table>	<input type="checkbox"/> Reagent Blank	<input type="checkbox"/> Auto Reagent Blank by Bottle	<input type="checkbox"/> Ask for calibration when reagent lot changes	<table><tr><td>R-Blank Limit</td></tr><tr><td>2500.0 mAbs</td></tr><tr><td>Cal Reps Range</td></tr><tr><td>15.0 %</td></tr><tr><td>Min Cal Reps</td></tr><tr><td>100.0 mAbs</td></tr><tr><td>Factor Change</td></tr><tr><td>35 %</td></tr><tr><td>M-Point Curve Fit</td></tr><tr><td>0.0 %</td></tr></table>	R-Blank Limit	2500.0 mAbs	Cal Reps Range	15.0 %	Min Cal Reps	100.0 mAbs	Factor Change	35 %	M-Point Curve Fit	0.0 %
Method	Curve Type																										
M-point	5P Log																										
Repeats	Stability																										
2																											
No																											
6																											
<input type="checkbox"/> Reagent Blank																											
<input type="checkbox"/> Auto Reagent Blank by Bottle																											
<input type="checkbox"/> Ask for calibration when reagent lot changes																											
R-Blank Limit																											
2500.0 mAbs																											
Cal Reps Range																											
15.0 %																											
Min Cal Reps																											
100.0 mAbs																											
Factor Change																											
35 %																											
M-Point Curve Fit																											
0.0 %																											

Calibrator		
	Calibrator	Conc
1	NaCl*	0.0
2	TC-CRP-1	**
3	TC-CRP-2	**
4	TC-CRP-3	**
5	TC-CRP-4	**
6	TC-CRP-5	**

# С-Реактивный белок, ЦРБ (CRP-FS)

ИттипoTD (Иммунотурбидиметрический)



\* -Рекомендуемый параметр.

\*\* - Вводится из паспорта калибратору.

Калибровка многоточечная (6 точек, 0+5 уровней) Spline или 5P Log,  
(калибратор TruCal CRP).

Контроль TruLab Protein уровень 1 и 2 или TruLab CRP уровень1 и 2