

Параметры для ввода в программу анализатора iLab-Taurus

Test Name: NN **Na** Test Code: IL_ **Na-d** Measure Sample Reagent Ranges Limits Calibration

Sample Type: **Serum**

Reaction Cycle: ☒ Standard ☐ Extended

Reporting Unit: ☒ ммоль/л * ☐ User Define
 Decimal Points: 1
 Conv.Factor: 0.000

Methodology:
 Type: ☐ End Point ☒ Rate
 Measuring Point: 22 - ☒ 29
 Photometric: ☐ 1 Wavelength ☒ 2 Wavelength
 Primary: 405 Secondary: 546
 Correction Constant: Slope: 1.000 Intercept: 0.000

Test Name: NN **Na** Test Code: IL_ **Na-d** Measure Sample Reagent Ranges Limits Calibration

Sample Volume:

	Sample	Dilution	
	Volume	Sample Vol.	Diluent Vol.
1	5.5	0.0	0.0
2	4.0	0.0	0.0
3	7.0	0.0	0.0
4			

Diluent:
 Diluent Warning Limit: 0 Tests

Reagent Volume:
 R1: **Na-d1** R2: **Na-d2**

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	120	0.0	<input checked="" type="checkbox"/>	20	No Control
R2	40	0.0	<input checked="" type="checkbox"/>	20	No Control

Sampling Condition:

	Condition No.	1	2	3	4
<input type="checkbox"/>	First Run	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Samp.Vol.Reduction	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Below N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Above N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Panic L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Panic H	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	User Range L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	User Range H	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Noise	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Prozone	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	HIGH!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	ABS!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

Ranges			Limits																																										
<table border="1"> <thead> <tr> <th colspan="3">Normal Range</th> </tr> <tr> <th></th> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>Male</td> <td>134</td> <td>150</td> </tr> <tr> <td>Female</td> <td>134</td> <td>150</td> </tr> <tr> <td>Other</td> <td>134</td> <td>150</td> </tr> </tbody> </table>			Normal Range				Lower	Upper	Male	134	150	Female	134	150	Other	134	150	<table border="1"> <thead> <tr> <th>Reaction Slope</th> <th>Absorbance Limit</th> </tr> </thead> <tbody> <tr> <td> <input type="radio"/> Negative <input checked="" type="radio"/> Positive </td> <td> <input checked="" type="radio"/> Above <input type="radio"/> Below </td> </tr> <tr> <td> <table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <td>25</td> <td>%</td> </tr> </thead></table> </td> <td> <table border="1"> <tr> <td>2500.0</td> <td>mAbs</td> </tr> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <td> <input type="radio"/> Above <input checked="" type="radio"/> Below </td> <td> <table border="1"> <tr> <td>Limit</td> <td>Equation</td> </tr> <tr> <td>0.0</td> <td>none</td> </tr> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <tr> <td>Judge Point</td> </tr> <tr> <td>0</td> </tr> </table> </td> </tr> </thead></table></td></tr></tbody> </table>		Reaction Slope	Absorbance Limit	<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below	<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <td>25</td> <td>%</td> </tr> </thead></table>	Non-Linear Limit		25	%	<table border="1"> <tr> <td>2500.0</td> <td>mAbs</td> </tr> </table>	2500.0	mAbs	<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <td> <input type="radio"/> Above <input checked="" type="radio"/> Below </td> <td> <table border="1"> <tr> <td>Limit</td> <td>Equation</td> </tr> <tr> <td>0.0</td> <td>none</td> </tr> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <tr> <td>Judge Point</td> </tr> <tr> <td>0</td> </tr> </table> </td> </tr> </thead></table>		Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below	<table border="1"> <tr> <td>Limit</td> <td>Equation</td> </tr> <tr> <td>0.0</td> <td>none</td> </tr> </table>	Limit	Equation	0.0	none	<table border="1"> <tr> <td>Judge Point</td> </tr> <tr> <td>0</td> </tr> </table>		Judge Point	0
Normal Range																																													
	Lower	Upper																																											
Male	134	150																																											
Female	134	150																																											
Other	134	150																																											
Reaction Slope	Absorbance Limit																																												
<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below																																												
<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <td>25</td> <td>%</td> </tr> </thead></table>	Non-Linear Limit		25	%	<table border="1"> <tr> <td>2500.0</td> <td>mAbs</td> </tr> </table>	2500.0	mAbs																																						
Non-Linear Limit																																													
25	%																																												
2500.0	mAbs																																												
<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <td> <input type="radio"/> Above <input checked="" type="radio"/> Below </td> <td> <table border="1"> <tr> <td>Limit</td> <td>Equation</td> </tr> <tr> <td>0.0</td> <td>none</td> </tr> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <tr> <td>Judge Point</td> </tr> <tr> <td>0</td> </tr> </table> </td> </tr> </thead></table>		Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below	<table border="1"> <tr> <td>Limit</td> <td>Equation</td> </tr> <tr> <td>0.0</td> <td>none</td> </tr> </table>	Limit	Equation	0.0	none	<table border="1"> <tr> <td>Judge Point</td> </tr> <tr> <td>0</td> </tr> </table>		Judge Point	0																																
Prozone Limit																																													
<input type="radio"/> Above <input checked="" type="radio"/> Below	<table border="1"> <tr> <td>Limit</td> <td>Equation</td> </tr> <tr> <td>0.0</td> <td>none</td> </tr> </table>	Limit	Equation	0.0	none																																								
Limit	Equation																																												
0.0	none																																												
<table border="1"> <tr> <td>Judge Point</td> </tr> <tr> <td>0</td> </tr> </table>		Judge Point	0																																										
Judge Point																																													
0																																													
<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>110</td> <td>180</td> </tr> </tbody> </table>			User Range		Lower	Upper	110	180																																					
User Range																																													
Lower	Upper																																												
110	180																																												
<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>110</td> <td>180</td> </tr> </tbody> </table>			ValidRange		Lower	Upper	110	180																																					
ValidRange																																													
Lower	Upper																																												
110	180																																												
<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>			Serum Index Limits			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0																																	
Serum Index Limits																																													
	Limit																																												
Hemolysis	0.0																																												
Icterus	0.0																																												
Lipemia	0.0																																												
<table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <td> <input type="radio"/> On <input checked="" type="radio"/> Off </td> <td></td> </tr> </thead></table>			Qualitative		<input type="radio"/> On <input checked="" type="radio"/> Off																																								
Qualitative																																													
<input type="radio"/> On <input checked="" type="radio"/> Off																																													

Test Name		Test Code		Measure	Sample Reagent	Ranges Limits	Calibration																															
NN	Na	IL	Na-d																																			
Calibration		Condition			Calibration/R-Blank Limit																																	
<table border="1"> <thead> <tr> <th>Method</th> <th>Curve Type</th> </tr> </thead> <tbody> <tr> <td>M-point</td> <td>Linear</td> </tr> <tr> <td>Repeats</td> <td>Stability</td> </tr> <tr> <td>2</td> <td></td> </tr> <tr> <td>No</td> <td></td> </tr> </tbody> </table>		Method	Curve Type	M-point	Linear	Repeats	Stability	2		No		<table border="1"> <tbody> <tr> <td><input type="checkbox"/></td> <td>Reagent Blank</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Auto Reagent Blank by Bottle</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Ask for calibration when reagent lot changes</td> </tr> </tbody> </table>			<input type="checkbox"/>	Reagent Blank	<input type="checkbox"/>	Auto Reagent Blank by Bottle	<input type="checkbox"/>	Ask for calibration when reagent lot changes	<table border="1"> <tbody> <tr> <td>R-Blank Limit</td> <td>2500.0</td> <td>mAbs</td> </tr> <tr> <td>Cal Reps Range</td> <td>20.0</td> <td>%</td> </tr> <tr> <td>Min Cal Reps</td> <td>100.0</td> <td>mAbs</td> </tr> <tr> <td>Factor Change</td> <td>20</td> <td>%</td> </tr> <tr> <td>M-Point Curve Fit</td> <td>0.0</td> <td>%</td> </tr> </tbody> </table>			R-Blank Limit	2500.0	mAbs	Cal Reps Range	20.0	%	Min Cal Reps	100.0	mAbs	Factor Change	20	%	M-Point Curve Fit	0.0	%
Method	Curve Type																																					
M-point	Linear																																					
Repeats	Stability																																					
2																																						
No																																						
<input type="checkbox"/>	Reagent Blank																																					
<input type="checkbox"/>	Auto Reagent Blank by Bottle																																					
<input type="checkbox"/>	Ask for calibration when reagent lot changes																																					
R-Blank Limit	2500.0	mAbs																																				
Cal Reps Range	20.0	%																																				
Min Cal Reps	100.0	mAbs																																				
Factor Change	20	%																																				
M-Point Curve Fit	0.0	%																																				
<table border="1"> <thead> <tr> <th colspan="3">Calibrator</th> </tr> <tr> <th></th> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Std-Na-1</td> <td>**</td> </tr> <tr> <td>2</td> <td>Std-Na-2</td> <td>**</td> </tr> <tr> <td>3</td> <td></td> <td></td> </tr> <tr> <td>4</td> <td></td> <td></td> </tr> <tr> <td>5</td> <td></td> <td></td> </tr> </tbody> </table>		Calibrator				Calibrator	Conc	1	Std-Na-1	**	2	Std-Na-2	**	3			4			5																		
Calibrator																																						
	Calibrator	Conc																																				
1	Std-Na-1	**																																				
2	Std-Na-2	**																																				
3																																						
4																																						
5																																						

Натрий

Ферментативный



* -Рекомендуемый параметр.

** - Вводится из паспорта калибратору.

Тип калибровки: Линейная, по двум стандартам из набора.

Контроль по TruLab N и TruLab P.