

Параметры для ввода в программу анализатора iLab-Taurus

Test Name: NN **C3c** Test Code: IL_ **C3c** Measure Sample Reagent Ranges Limits Calibration

Sample Type
Serum

Reaction Cycle
☒ Standard ☐ Extended

Reporting Unit
☒ г/л * ☐ User Define Decimal Points: 3
Conv. Factor: 0.000

Methodology
Type
☒ End Point ☐ Rate
Measuring Point
16 - ☒ 33
Photometric
☐ 1 Wavelength ☒ 2 Wavelength
Primary: 340 **Secondary**: 700
Correction Constant
Slope: 1.000 **Intercept**: 0.000

Test Name: NN **C3c** Test Code: IL_ **C3c** Measure Sample Reagent Ranges Limits Calibration

Sample Volume

	Sample Volume	Dilution Sample Vol.	Dilution Diluent Vol.
1	2.25	0.0	0.0
2	2.0	30.0	60.0
3	4.0	0.0	0.0
4	0.0	0.0	0.0

Diluent: **Diluent Warning Limit**: 0 Tests

Reagent Volume
R1: C3c-1 R2: C3c-2

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	150	0.0	<input checked="" type="checkbox"/>	20	No Control
R2	30	0.0	<input checked="" type="checkbox"/>	20	No Control

Sampling Condition

Condition No.	1	2	3	4
First Run	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Samp.Vol.Reduction	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Below N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Above N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Panic L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Panic H	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
User Range L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/> User Range H	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Noise	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Prozone	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/> HIGH!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/> ABS!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

Ranges		Limits																																																																																						
<table border="1"> <thead> <tr> <th colspan="3">Normal Range</th> </tr> <tr> <th></th> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>Male</td> <td>0.9</td> <td>1.8</td> </tr> <tr> <td>Female</td> <td>0.9</td> <td>1.8</td> </tr> <tr> <td>Other</td> <td>0.9</td> <td>1.8</td> </tr> </tbody> </table>		Normal Range				Lower	Upper	Male	0.9	1.8	Female	0.9	1.8	Other	0.9	1.8	<table border="1"> <thead> <tr> <th>Reaction Slope</th> <th>Absorbance Limit</th> </tr> </thead> <tbody> <tr> <td> <input type="radio"/> Negative <input checked="" type="radio"/> Positive </td> <td> <input checked="" type="radio"/> Above <input type="radio"/> Below </td> </tr> <tr> <td> <table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>100</td> <td></td> </tr> </tbody> </table> </td> <td> <table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2600.0</td> <td></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>5.00</td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th>Index</th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>5.00</td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <th>On</th> <th>Off</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/></td> <td><input checked="" type="radio"/></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table> </td> </tr> </tbody> </table>		Reaction Slope	Absorbance Limit	<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below	<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>100</td> <td></td> </tr> </tbody> </table>	Non-Linear Limit		Value	%	100		<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2600.0</td> <td></td> </tr> </tbody> </table>	Prozone Limit		Value	%	2600.0		<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>5.00</td> </tr> </tbody> </table>		User Range		Lower	Upper	0.01	5.00	<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th>Index</th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>		Serum Index Limits		Index	Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0	<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>5.00</td> </tr> </tbody> </table>		ValidRange		Lower	Upper	0.01	5.00	<table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table>		Equation		Limit	Equation	0.0	none	<table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <th>On</th> <th>Off</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/></td> <td><input checked="" type="radio"/></td> </tr> </tbody> </table>		Qualitative		On	Off	<input type="radio"/>	<input checked="" type="radio"/>	<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table>		Judge Point		Value	%	0	
Normal Range																																																																																								
	Lower	Upper																																																																																						
Male	0.9	1.8																																																																																						
Female	0.9	1.8																																																																																						
Other	0.9	1.8																																																																																						
Reaction Slope	Absorbance Limit																																																																																							
<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below																																																																																							
<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>100</td> <td></td> </tr> </tbody> </table>	Non-Linear Limit		Value	%	100		<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2600.0</td> <td></td> </tr> </tbody> </table>	Prozone Limit		Value	%	2600.0																																																																												
Non-Linear Limit																																																																																								
Value	%																																																																																							
100																																																																																								
Prozone Limit																																																																																								
Value	%																																																																																							
2600.0																																																																																								
<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>5.00</td> </tr> </tbody> </table>		User Range		Lower	Upper	0.01	5.00	<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th>Index</th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>		Serum Index Limits		Index	Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0																																																																					
User Range																																																																																								
Lower	Upper																																																																																							
0.01	5.00																																																																																							
Serum Index Limits																																																																																								
Index	Limit																																																																																							
Hemolysis	0.0																																																																																							
Icterus	0.0																																																																																							
Lipemia	0.0																																																																																							
<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>5.00</td> </tr> </tbody> </table>		ValidRange		Lower	Upper	0.01	5.00	<table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <th>Limit</th> <th>Equation</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td>none</td> </tr> </tbody> </table>		Equation		Limit	Equation	0.0	none																																																																									
ValidRange																																																																																								
Lower	Upper																																																																																							
0.01	5.00																																																																																							
Equation																																																																																								
Limit	Equation																																																																																							
0.0	none																																																																																							
<table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> <tr> <th>On</th> <th>Off</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/></td> <td><input checked="" type="radio"/></td> </tr> </tbody> </table>		Qualitative		On	Off	<input type="radio"/>	<input checked="" type="radio"/>	<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td></td> </tr> </tbody> </table>		Judge Point		Value	%	0																																																																										
Qualitative																																																																																								
On	Off																																																																																							
<input type="radio"/>	<input checked="" type="radio"/>																																																																																							
Judge Point																																																																																								
Value	%																																																																																							
0																																																																																								

Test Name		Test Code		Measure		Sample Reagent		Ranges Limits		Calibration																																							
NN	C3c	IL	C3c																																														
Calibration				Condition				Calibration/R-Blank Limit																																									
<table border="1"> <thead> <tr> <th colspan="2">Method</th> <th colspan="2">Curve Type</th> </tr> </thead> <tbody> <tr> <td>M-point</td> <td></td> <td>5P Log</td> <td></td> </tr> </tbody> </table>				Method		Curve Type		M-point		5P Log		<table border="1"> <thead> <tr> <th colspan="2">Repeats</th> <th colspan="2">Stability</th> </tr> </thead> <tbody> <tr> <td>2</td> <td></td> <td></td> <td>days</td> </tr> </tbody> </table>				Repeats		Stability		2			days	<table border="1"> <thead> <tr> <th colspan="2">No</th> </tr> </thead> <tbody> <tr> <td>6</td> <td></td> </tr> </tbody> </table>				No		6																			
Method		Curve Type																																															
M-point		5P Log																																															
Repeats		Stability																																															
2			days																																														
No																																																	
6																																																	
<table border="1"> <thead> <tr> <th colspan="3">Calibrator</th> </tr> <tr> <th></th> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>NaCl*</td> <td>0.0</td> </tr> <tr> <td>2</td> <td>TC-prot-1</td> <td>**</td> </tr> <tr> <td>3</td> <td>TC-prot-2</td> <td>**</td> </tr> <tr> <td>4</td> <td>TC-prot-3</td> <td>**</td> </tr> <tr> <td>5</td> <td>TC-prot-4</td> <td>**</td> </tr> <tr> <td>6</td> <td>TC-prot-5</td> <td>**</td> </tr> </tbody> </table>				Calibrator				Calibrator	Conc	1	NaCl*	0.0	2	TC-prot-1	**	3	TC-prot-2	**	4	TC-prot-3	**	5	TC-prot-4	**	6	TC-prot-5	**	<table border="1"> <thead> <tr> <th colspan="2">Condition</th> </tr> </thead> <tbody> <tr> <td><input type="checkbox"/></td> <td>Reagent Blank</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Auto Reagent Blank by Bottle</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Ask for calibration when reagent lot changes</td> </tr> </tbody> </table>				Condition		<input type="checkbox"/>	Reagent Blank	<input type="checkbox"/>	Auto Reagent Blank by Bottle	<input type="checkbox"/>	Ask for calibration when reagent lot changes	<table border="1"> <thead> <tr> <th colspan="2">R-Blank Limit</th> </tr> <tr> <th>Value</th> <th>mAbs</th> </tr> </thead> <tbody> <tr> <td>2500.0</td> <td></td> </tr> </tbody> </table>				R-Blank Limit		Value	mAbs	2500.0	
Calibrator																																																	
	Calibrator	Conc																																															
1	NaCl*	0.0																																															
2	TC-prot-1	**																																															
3	TC-prot-2	**																																															
4	TC-prot-3	**																																															
5	TC-prot-4	**																																															
6	TC-prot-5	**																																															
Condition																																																	
<input type="checkbox"/>	Reagent Blank																																																
<input type="checkbox"/>	Auto Reagent Blank by Bottle																																																
<input type="checkbox"/>	Ask for calibration when reagent lot changes																																																
R-Blank Limit																																																	
Value	mAbs																																																
2500.0																																																	
				<table border="1"> <thead> <tr> <th colspan="2">Cal Reps Range</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>15.0</td> <td></td> </tr> </tbody> </table>				Cal Reps Range		Value	%	15.0																																					
Cal Reps Range																																																	
Value	%																																																
15.0																																																	
				<table border="1"> <thead> <tr> <th colspan="2">Min Cal Reps</th> </tr> <tr> <th>Value</th> <th>mAbs</th> </tr> </thead> <tbody> <tr> <td>200.0</td> <td></td> </tr> </tbody> </table>				Min Cal Reps		Value	mAbs	200.0																																					
Min Cal Reps																																																	
Value	mAbs																																																
200.0																																																	
				<table border="1"> <thead> <tr> <th colspan="2">Factor Change</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>35</td> <td></td> </tr> </tbody> </table>				Factor Change		Value	%	35																																					
Factor Change																																																	
Value	%																																																
35																																																	
				<table border="1"> <thead> <tr> <th colspan="2">M-Point Curve Fit</th> </tr> <tr> <th>Value</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0.0</td> <td></td> </tr> </tbody> </table>				M-Point Curve Fit		Value	%	0.0																																					
M-Point Curve Fit																																																	
Value	%																																																
0.0																																																	

* -Рекомендуемый параметр.

** - Вводится из паспорта калибратору.

Калибровка многоточечная (6 точек, 0+5 уровней) Spline или 5P Log, калибратор TruCal Protein.

Контроль TruLab Protein уровень 1 и 2.