

# Щелочная фосфатаза

DGKC-1970 (п-нитрофенольный метод)

## Параметры для ввода в программу анализатора iLab-Taurus

Test Name: NN ALP Test Code: IL\_ ALP-d

Measure Sample Reagent Ranges Limits Calibration

Sample Type: Serum

Reaction Cycle: Standard Extended

Reporting Unit: Ед/л \* Decimal Points: 0

User Define: Conv.Factor: 0.000

Methodology: Type: End Point Rate

Measuring Point: 19 - 29

Photometric: 1 Wavelength 2 Wavelength

Primary: 405 Secondary: 546

Correction Constant: Slope: 1.000 Intercept: 0.000

Test Name: NN ALP Test Code: IL\_ ALP-d

Measure Sample Reagent Ranges Limits Calibration

Sample Volume: Sampling Condition:

Sample	Dilution	Diluent	Condition No.	1	2	3	4
1	3.0	0.0	First Run				
2	3.0	10.0	Samp.Vol.Reduction				
3	6.0	0.0	Below N-Range				
4	2.0	5.0	Above N-Range				
			Panic L				
			Panic H				
			User Range L				
			User Range H				
			Noise				
			Prozone				
			HIGH!				
			ABS!				

Diluent: Diluent Warning Limit: 0 Tests

Reagent Volume: R1: aAmy-d1 R2: aAmy-d2

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	140	0.0		20	No Control
R2	35	0.0		20	No Control

# Щелочная фосфатаза

DGKC-1970 (п-нитрофенольный метод)

Ranges		Limits																																												
<table border="1"><thead><tr><th colspan="3">Normal Range</th></tr><tr><th></th><th>Lower</th><th>Upper</th></tr></thead><tbody><tr><td>Male</td><td>0</td><td>258</td></tr><tr><td>Female</td><td>0</td><td>258</td></tr><tr><td>Other</td><td>0</td><td>258</td></tr></tbody></table>		Normal Range				Lower	Upper	Male	0	258	Female	0	258	Other	0	258	<table border="1"><thead><tr><th>Reaction Slope</th><th>Absorbance Limit</th></tr></thead><tbody><tr><td><input type="radio"/> Negative <input checked="" type="radio"/> Positive</td><td><input checked="" type="radio"/> Above <input type="radio"/> Below</td></tr><tr><td><table border="1"><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>25</td><td>%</td></tr></thead></table></td><td><table border="1"><thead><tr><td>2500.0</td><td>mAbs</td></tr></thead></table></td></tr><tr><td colspan="2"><table border="1"><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td><table border="1"><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table></td><td></td></tr><tr><td colspan="2"><table border="1"><thead><tr><th>Judge Point</th></tr><tr><td>0</td></tr></thead></table></td></tr></thead></table></td></tr></tbody></table>		Reaction Slope	Absorbance Limit	<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below	<table border="1"><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>25</td><td>%</td></tr></thead></table>	Non-Linear Limit		25	%	<table border="1"><thead><tr><td>2500.0</td><td>mAbs</td></tr></thead></table>	2500.0	mAbs	<table border="1"><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td><table border="1"><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table></td><td></td></tr><tr><td colspan="2"><table border="1"><thead><tr><th>Judge Point</th></tr><tr><td>0</td></tr></thead></table></td></tr></thead></table>		Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below		<table border="1"><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table>	Limit	Equation	0.0	none		<table border="1"><thead><tr><th>Judge Point</th></tr><tr><td>0</td></tr></thead></table>		Judge Point	0
Normal Range																																														
	Lower	Upper																																												
Male	0	258																																												
Female	0	258																																												
Other	0	258																																												
Reaction Slope	Absorbance Limit																																													
<input type="radio"/> Negative <input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above <input type="radio"/> Below																																													
<table border="1"><thead><tr><th colspan="2">Non-Linear Limit</th></tr><tr><td>25</td><td>%</td></tr></thead></table>	Non-Linear Limit		25	%	<table border="1"><thead><tr><td>2500.0</td><td>mAbs</td></tr></thead></table>	2500.0	mAbs																																							
Non-Linear Limit																																														
25	%																																													
2500.0	mAbs																																													
<table border="1"><thead><tr><th colspan="2">Prozone Limit</th></tr><tr><td><input type="radio"/> Above <input checked="" type="radio"/> Below</td><td></td></tr><tr><td><table border="1"><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table></td><td></td></tr><tr><td colspan="2"><table border="1"><thead><tr><th>Judge Point</th></tr><tr><td>0</td></tr></thead></table></td></tr></thead></table>		Prozone Limit		<input type="radio"/> Above <input checked="" type="radio"/> Below		<table border="1"><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table>	Limit	Equation	0.0	none		<table border="1"><thead><tr><th>Judge Point</th></tr><tr><td>0</td></tr></thead></table>		Judge Point	0																															
Prozone Limit																																														
<input type="radio"/> Above <input checked="" type="radio"/> Below																																														
<table border="1"><thead><tr><th>Limit</th><th>Equation</th></tr></thead><tbody><tr><td>0.0</td><td>none</td></tr></tbody></table>	Limit	Equation	0.0	none																																										
Limit	Equation																																													
0.0	none																																													
<table border="1"><thead><tr><th>Judge Point</th></tr><tr><td>0</td></tr></thead></table>		Judge Point	0																																											
Judge Point																																														
0																																														

User Range	
Lower	Upper
3	4500

ValidRange	
Lower	Upper
3	4500

Qualitative	
<input type="radio"/> On <input checked="" type="radio"/> Off	

Serum Index Limits	
	Limit
Hemolysis	0.0
Icterus	0.0
Lipemia	0.0

Test Name		Test Code		Measure		Sample Reagent		Ranges Limits		Calibration																																																							
NN	ALP	IL_	ALP-d																																																														
Calibration				Condition				Calibration/R-Blank Limit																																																									
<table border="1"><thead><tr><th colspan="2">Method</th></tr><tr><td>1-point</td><td></td></tr></thead><tbody><tr><td colspan="2">Curve Type</td></tr><tr><td>Linear</td><td></td></tr><tr><td colspan="2">Repeats</td></tr><tr><td>2</td><td></td></tr><tr><td colspan="2">Stability</td></tr><tr><td></td><td>days</td></tr><tr><td colspan="2">No</td></tr><tr><td></td><td></td></tr></tbody></table>				Method		1-point		Curve Type		Linear		Repeats		2		Stability			days	No				<table border="1"><thead><tr><th colspan="2">Calibrator</th></tr><tr><th>Calibrator</th><th>Conc</th></tr></thead><tbody><tr><td>1</td><td>TruCalU **</td></tr><tr><td>2</td><td></td></tr><tr><td>3</td><td></td></tr><tr><td>4</td><td></td></tr><tr><td>5</td><td></td></tr></tbody></table>				Calibrator		Calibrator	Conc	1	TruCalU **	2		3		4		5		<table border="1"><thead><tr><th colspan="2">R-Blank Limit</th></tr><tr><td>2500.0</td><td>mAbs</td></tr><tr><td colspan="2">Cal Reps Range</td></tr><tr><td>20.0</td><td>%</td></tr><tr><td colspan="2">Min Cal Reps</td></tr><tr><td>50.0</td><td>mAbs</td></tr><tr><td colspan="2">Factor Change</td></tr><tr><td>20</td><td>%</td></tr><tr><td colspan="2">M-Point Curve Fit</td></tr><tr><td>0.0</td><td>%</td></tr></thead></table>				R-Blank Limit		2500.0	mAbs	Cal Reps Range		20.0	%	Min Cal Reps		50.0	mAbs	Factor Change		20	%	M-Point Curve Fit		0.0	%
Method																																																																	
1-point																																																																	
Curve Type																																																																	
Linear																																																																	
Repeats																																																																	
2																																																																	
Stability																																																																	
	days																																																																
No																																																																	
Calibrator																																																																	
Calibrator	Conc																																																																
1	TruCalU **																																																																
2																																																																	
3																																																																	
4																																																																	
5																																																																	
R-Blank Limit																																																																	
2500.0	mAbs																																																																
Cal Reps Range																																																																	
20.0	%																																																																
Min Cal Reps																																																																	
50.0	mAbs																																																																
Factor Change																																																																	
20	%																																																																
M-Point Curve Fit																																																																	
0.0	%																																																																

| | Condition                           |  | |-------------------------------------|--| | <input checked="" type="checkbox"/> | Reagent Blank                                | | <input type="checkbox"/>            | Auto Reagent Blank by Bottle                 | | <input type="checkbox"/>            | Ask for calibration when reagent lot changes | | | | |  | | | |

# Щелочная фосфатаза

DGKC-1970 (*n*-нитрофенольный метод)

---



\* -Рекомендуемый параметр.

\*\* - Вводится из паспорта калибратору.

Тип калибровки: Линейная (калибратор TruCalU) или по фактору (см инструкцию).  
Контроль по TruLab N и TruLab P.