

Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

☒ General
 ☐ Calibration
 ☐ SmartWash
 ☐ Results
 ☐ Interpretation

Assay:
 Type: Photometric
 Version: 1

Number:
 Assay availability:
 Date: 18.02.2020

Run controls for onboard reagents by:
 Time: 11:01:08

Operator: ADMIN

☒ Reaction definition
 ☐ Reagent / Sample
 ☐ Validity checks

Reaction mode:

Primary Secondary Read times

Wavelength: /
 Main: --

Last required read:
 Flex: --

Absorbance range: --
 Color correction: --

Sample blank type:

Configure assay parameters...

☒ General
 ☐ Calibration
 ☐ SmartWash
 ☐ Results
 ☐ Interpretation

Assay:
 Type: Photometric
 Version: 1

Number:
 Assay availability:
 Date: 20.02.2020

Run controls for onboard reagents by:
 Time: 17:14:07

Operator: FSE

☐ Reaction definition
 ☒ Reagent / Sample
 ☐ Validity checks

Reagent:
 R1 R2

Reagent volume:

Diluent:
 Water volume:

Diluent dispense mode:
 Dispense mode:

| Dilution name | Sample | Diluted sample | Diluent | Water | Dilution factor | Default dilution |
|---------------|--------|----------------|---------|-------|-----------------|----------------------------------|
| Standart | 4.0 | | | | 1:1.00 | <input checked="" type="radio"/> |
| 1:5 | 20.0 | 4.0 | 80.0 | | 1:5.00 | <input type="radio"/> |
| 2x | 8.0 | | | | 2:1.00 | <input type="radio"/> |

(Myogl)

ИммуноTD (Иммунотурбидиметрический, усиленный латексом)

Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Myogl Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: Calibrator level: Concentration:

TruCalMyo Blank: Water 0.0000

Cal 1: TruCalMyo1
Cal 2: TruCalMyo2
Cal 3: TruCalMyo3
Cal 4: TruCalMyo4
Cal 5: None
Cal 6: None

Replicates: 2 [Range 1 - 3]

Assay of

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Myogl Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator:

| Calibrator level | Sample | Diluted sample | Diluent | Water |
|-------------------|--------|----------------|---------|-------|
| Blank: Water | 4.0 | | | |
| Cal 1: TruCalMyo1 | 4.0 | | | |
| Cal 2: TruCalMyo2 | 4.0 | | | |
| Cal 3: TruCalMyo3 | 4.0 | | | |
| Cal 4: TruCalMyo4 | 4.0 | | | |
| Cal 5: None | | | | |
| Cal 6: None | | | | |

Assay of

Done Cancel

?

Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

☐ General ☒ Calibration ☐ SmartWash ☐ Results ☐ Interpretation

Assay: Myogl Assay number: *** Date: 20.02.2020
Calibration method: Spline Time: 17:14
Operator: FSE

☐ Calibrators ☐ Volumes ☒ Intervals ☐ Validity checks

Calibration intervals:

Full interval: 0 (hours)
Adjust interval: 0 (hours)

Calibration type:

Adjust type: 2 Point
Adjust level: Cal 1
Default ordering type: Full

Assay of

Done
Cancel

?

Единицы измерения - мкг/л; знаков после запятой 1 (**Configure result units**).

Диапазон метода: 1-600 мкг/л (**Configure assay parameters\ Results**).

Референсный интервал: 0-70 мкг/л (**Configure assay parameters\ Results**).

При наличии ЛИС не обязательно.

Диапазон нормальных значений указан в соответствии с рекомендациями производителя реагентов.

Калибратор TruCal Myoglobin, набор из 4-х уровней, нулевой - вода (или физраствор).

Контроль по TruLab Protein уровень 1 и уровень 2.