

## Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay:  Type: Photometric Version: 1  
 Number:  Assay availability:  Date: 18.02.2020  
 Run controls for onboard reagents by:  Time: 11:01:08  
 Operator: ADMIN

Reaction definition Reagent / Sample Validity checks

Reaction mode:

Primary Secondary Read times  
 Wavelength:  /  Main:  --   
 Last required read:   
 Absorbance range:  --  Color correction:  --   
 Sample blank type:  Blank:  --

Assay 13 of 54

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay:  Type: Photometric Version: 1  
 Number:  Assay availability:  Date: 20.02.2020  
 Run controls for onboard reagents by:  Time: 17:14:07  
 Operator: FSE

Reaction definition Reagent / Sample Validity checks

Reagent:   
 Diluent:   
 Diluent dispense mode:   
 Reagent volume:  R1  
 Water volume:  R2  
 Dispense mode:

Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution
Standart	12.0				1:1.00	<input checked="" type="radio"/>
1:3	4.0				1:2.89	<input type="radio"/>
						<input type="radio"/>

Assay of

Done Cancel

?

## Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: UBIC Assay number: \*\*\* Date: 20.02.2020  
 Calibration method: Linear Time: 17:14  
 Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: TruCal\_U  
 Replicates: 2 [Range 1 - 3]

Calibrator level: Blank: Water  
 Cal 1: TruCal\_U1  
 Cal 2: None  
 Cal 3: None  
 Cal 4: None  
 Cal 5: None  
 Cal 6: None

Concentration: 0.0000

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: UBIC Assay number: \*\*\* Date: 20.02.2020  
 Calibration method: Linear Time: 17:14  
 Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator:

Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	12.0			
Cal 1: TruCal_U1	12.0			
Cal 2: None				
Cal 3: None				
Cal 4: None				
Cal 5: None				
Cal 6: None				

Done Cancel

?

---

Параметры для ввода в программу анализатора Architect c8000 / c4000

Единицы измерения - мкмоль /л; знаков после запятой 1 (**Configure result units**).

Диапазон метода: 1-135 мкмоль/л (**Configure assay parameters\ Results**).

Референсный интервал: 21-84 мкмоль/л (**Configure assay parameters\ Results**).

При наличии ЛИС, необязательно.

Диапазон нормальных значений указан в соответствии с рекомендациями производителя реагентов.

Калибровка линейная, калибратор TruCal U, нулевой - вода (или физраствор).

Контроль по TruLab N.