

Комплемент C4

ИммуноTD (Иммунотурбидиметрический)



Параметры для ввода в программу анализатора Beckman-Coulter AU 480

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <input type="text" value="C4"/> < > Type: <input type="text" value="Serum"/> Operation <input type="text" value="Yes"/>					
Sample Volume <input type="text" value="2.6"/> μL		Dilution <input type="text" value="0"/> μL		OD Limit	
Pre-Dilution Rate <input type="text" value="1"/>		Dilution <input type="text" value="0"/> μL		Min.OD <input type="text" value="-2.0"/> Max.OD <input type="text" value="2.5"/>	
Rgt. Volume R1(R1-1) <input type="text" value="90"/> μL		Dilution <input type="text" value="0"/> μL		Reagent OD Limit	
				First Low <input type="text" value="-2.0"/> High <input type="text" value="2.5"/>	
				Last Low <input type="text" value="-2.0"/> High <input type="text" value="2.5"/>	
R2(R2-1) <input type="text" value="18"/> μL		Dilution <input type="text" value="0"/> μL		Dynamic Range Low <input type="text" value="0"/> High <input type="text" value="0.9"/>	
				Correlation Factor A <input type="text" value="1"/> B <input type="text" value="0"/>	
Wavelength Pri <input type="text" value="340"/> nm		Sec. <input type="text" value="800"/> nm		Factor for Maker A <input type="text" value="1"/> B <input type="text" value="0"/>	
Method <input type="text" value="END"/>		Onboard Stability Period <input type="text" value="999"/> Day <input type="text" value=""/> Hour			
Reaction Slope <input type="text" value="+"/> ∇					
Measuring Point1 First <input type="text" value="0"/>		Last <input type="text" value="27"/>			
Measuring Point2 First <input type="text" value="0"/>		Last <input type="text" value="10"/>			
Linearity Limit <input type="text" value=""/> %					
Lag Time Check <input type="text" value=""/>					

Parameters		Specific Test Parameters			
General	LIH	ISE	HbA1c	Calculated Test	Range
Test Name: <input type="text" value="C4"/> < > Type: <input type="text" value="Serum"/>					
Value/Flag: <input type="text" value=""/>		Level L: <input type="text" value=""/>		Level H: <input type="text" value=""/>	
Specific Ranges:					
		From To		Low High	
Sex Year Month Year Month					
1. <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/>					
2. <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/>					
3. <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/>					
4. <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/>					
5. <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/>					
6. <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/> <input type="text" value=""/>					
7. No demographics				<input type="text" value="0.1"/> <input type="text" value="0.4"/>	
8. Not within expected values				<input type="text" value="0.1"/> <input type="text" value="0.4"/>	
Unit <input type="text" value="g/L"/>		Decimal Places <input type="text" value="2"/>			
Panic Value					
Low <input type="text" value=""/> High <input type="text" value=""/>					

Parameters		Calibration Parameters	
Calibrators	Calibration Specific	STAT Table Calibration	
General	ISE		
Test Name: <input type="text" value="C4"/> < > Type: <input type="text" value="Serum"/> <input type="radio"/> Use Serum Cal.			
Calibration Type: <input type="text" value="6AB"/>		Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="2"/>	
<Calibrator Parameters>			
Calibrator	OD	Conc	Factor/OD range Low High
Point 1: TruCal Prot-0		*	-0.1 2.5
Point 2: TruCal Prot-1		*	-0.1 2.5
Point 3: TruCal Prot-2		*	-0.1 2.5
Point 4: TruCal Prot-3		*	-0.1 2.5
Point 5: TruCal Prot-4		*	-0.1 2.5
Point 6: TruCal Prot-5		*	-0.1 2.5
Point 7:			
Point 8:			
Point 9:			
Point 10:			
Slope Check <input type="text" value="None"/>		Allowance Range Check	
<input type="radio"/> Reagent Blank <input type="text" value=""/>		<input type="radio"/> Calibration <input type="text" value=""/>	
Advanced Calibration Operation <input type="text" value=""/>		Interval (RB/ACAL) <input type="text" value=""/>	
<Point Cal. For No. of Correction Points <input type="text" value=""/>		Use Master Curve <input type="text" value=""/>	
Master Curve		<input type="radio"/> Lot Calibration	
Calibrator	OD	Conc	Low High
Point 1:			
Point 2:			
Stability		Reagent Blank <input type="text" value=""/> Day <input type="text" value=""/> Hour	
		Calibration <input type="text" value=""/> Day <input type="text" value=""/> Hour	
MB Type Factor: <input type="text" value=""/>		1-Point Calibration Point <input type="text" value=""/>	
		<input type="radio"/> with Conc-0	

Диапазон нормальных значений указан в соответствии с рекомендациями производителя.

При использовании единиц измерения, отличающихся от приведенных, убедитесь, что значения стандартов, контрольных материалов, диапазонов нормальных значений и линейности метода введены в этих же единицах.

*-вводится в соответствии с настройками системы

** -вводится из паспорта к калибратору (TruCal Protein).

Контроль по TruLab Protein уровень 1 и уровень 2.