

TruCal Lipid

LOT 35402

Cat. No. 1 3570 99 10 45 (3 x 2 mL)

ВНИМАНИЕ!
Уточняйте паспортные значения
на www.diakonlab.ru или
по тел.: 8-800-200-63-39



2026-08-31



Updated sheets: www.diasys-diagnostics.com/service-area/value-sheets/

Value sheet version: 2

Short Name / Analyte	Method	Reagent ID	Traceability	Value	Unit
HDL-C HDL cholesterol Холестерин-ЛПВП	Immunoinhibition	1 3521	NIST SRM 1951 Level 2	64,5	mg/dL
	Прямой иммунологический			1,67	mmol/L
				0,645	g/L
	Direct, colorimetric	1 3561	Commercially available	60,9	mg/dL
	Прямой колориметрический		assay traceable to CDC	1,57	mmol/L
			reference method	0,609	g/L
LDL-C LDL cholesterol Холестерин-ЛПНП	Select	1 4121	NIST SRM 1951 Level 2	140	mg/dL
	Прямой селективный			3,61	mmol/L
				1,40	g/L
	Direct, colorimetric	1 4131	NIST SRM 1951 Level 2	149	mg/dL
	Прямой колориметрический			3,85	mmol/L
				1,49	g/L
NEFA Non esterified fatty acids Свободные жирные кислоты	Enzymatic	1 5781	Internal Standard	51,6	mg/dL
	Ферментативный			1,83	mmol/L
				0,516	g/L
PL Phospholipids Фосфолипиды	Enzymatic	1 5741	Internal Standard	318	mg/dL
	Ферментативный			4,10	mmol/L
				3,18	g/L



DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim Germany
www.diasys-diagnostics.com