


TruCal Lipid

LOT 34675

Cat. No. 1 3570 99 10 45 (3 x 2 mL)

ВНИМАНИЕ!
Уточняйте паспортные значения
на www.diakonlab.ru или
по тел.: 8-800-200-63-39

 **2026-04-30**

 **IVD**

Updated sheets: www.diasys-diagnostics.com/service-area/value-sheets/

Value sheet version: 2

Short Name / Analyte	Method	Reagent ID	Traceability	Value	Unit
HDL-C HDL cholesterol Холестерин-ЛПВП	Immunoinhibition	1 3521	NIST SRM 1951 Level 2	57,5	mg/dL
	Прямой иммунологический			1,49	mmol/L
				0,575	g/L
	Direct, colorimetric	1 3561	Commercially available	55,6	mg/dL
	Прямой колориметрический		assay traceable to CDC	1,44	mmol/L
			reference method	0,556	g/L
LDL-C LDL cholesterol Холестерин-ЛПНП	Select	1 4121	NIST SRM 1951 Level 2	140	mg/dL
	Прямой селективный			3,62	mmol/L
				1,40	g/L
	Direct, colorimetric	1 4131	NIST SRM 1951 Level 2	153	mg/dL
	Прямой колориметрический			3,96	mmol/L
				1,53	g/L
NEFA Non esterified fatty acids Свободные жирные кислоты	Enzymatic	1 5781	Internal Standard	52,5	mg/dL
	Ферментативный			1,86	mmol/L
				0,525	g/L
PL Phospholipids Фосфолипиды	Enzymatic	1 5741	Internal Standard	287	mg/dL
	Ферментативный			3,70	mmol/L
				2,87	g/L



DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim Germany
www.diasys-diagnostics.com